JXA-8230 JXA-8530F

SUPERPROBE

***EPMA
JXA-8230***

HYPERPROBE

***FE-EPMA
JXA-8530F***

JEOL

Serving Advanced Technology

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— JXA-8230/JXA-8530F —

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The screenshot displays a software interface for managing analysis recipes. It is divided into several sections:

- 面分析 (Surface Analysis):** Includes fields for recipe name (るつぼ), comments, and an '分析開始' (Start Analysis) button.
- 電子光学条件 (Electron Optics Conditions):** Shows parameters like '加速電圧[kV]' (15.0) and '照射電流[A]' (1.00 x 10⁻⁸), with a '読み込み' (Load) button.
- 位置条件 (Position Conditions):** A table listing analysis positions. A callout bubble points to the 'No. 2' entry with the text: **新たな分析位置を登録。**
- 拡張条件 (Expansion Conditions):** Includes '通常エリア指定' (Normal Area Designation) and '拡大率' (Magnification) set to 256 x 256.
- 表示項目 (Display Items):** A table showing analysis channels and their results.
- EDS/IMS Parameters:** Checkboxes for 'SEI', 'COMPO', 'TOPO', 'AUX1', 'AUX2', and 'AUX3'. The 'IMS' section is active, showing '動作モード' (Operation Mode) set to '現在値を使用' (Use Current Value) and 'コントラスト' (Contrast) and 'ブライツネス' (Brightness) both set to 0.

Additional callouts include: **元素条件が読み込まれる。 ※編集も可能** (Element conditions are loaded. ※Editing is also possible) pointing to the '読み込み' button, and **拡大率** (Magnification) pointing to the '256 x 256' field.

No.	コメント	Scan	PX	PY	時間[ms]	Size
0001	るつぼ	Beam	256	256	10.0	-
0002						

CH-1	CH-2	CH-3	CH-4	CH-5
Zr(PETH)	Al(TAP)	Sn(PETJ)	Fe(LIFH)	

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