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
Figure 1. Schematic diagram of the AccuTOF atmospheric pressure interface (API)


Figure 2. The white DART ceramic insulator cap is positioned approximately 1 cm from the apex of orifice 1 (the silver cone on the right).

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

Figure 3. Schematic diagram of VAPUR® interface. On the AccuTOF, the gap between orifice 1 and the exit of the VAPUR® ceramic tube should be 2 mm for optimal performance.

Figure 4. The VAPUR® interface mounted on the AccuTOF DART.

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
Figure 5. (a) positive-ion low-mass DART background without the VAPUR® installed, showing the dominant reagent ions with trace laboratory solvent peaks and (b) the low-mass background observed with the VAPUR® installed.

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
Figure 6. Comparison of signal for quinine, epitestosterone, and methyl stearate (a) without the VAPUR® interface and (b) with the VAPUR® interface.

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
Figure 7. Positive-ion DART mass spectra of a PEG sample measured (a) with no VAPUR® and (b) with the VAPUR® installed.

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
Figure 8. Chromatograms for three replicate measurements of a sample containing diisobutyl phthalate (DIBP) (a) with no VAPUR® and (b) with the VAPUR® installed.

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Copyright(C)2016 JEOL Ltd. All Rights Reserved.

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1. Cody, R.B., J.A. Laramée, and H.D. Durst, *Versatile New Ion Source for the Analysis of Materials in Open Air under Ambient Conditions*. Analytical Chemistry, 2005. **77**(8): p. 2297-2302.
2. Cody, R.B., et al., *Direct Analysis in Real Time (DART™) Mass Spectrometry*. JEOL News, 2005. **40**(1): p. 8-12.
3. Fernández, F.M., et al., *Characterization of Solid Counterfeit Drug Samples by Desorption Electrospray Ionization and Direct-analysis-in-real-time Coupled to Time-of-flight Mass Spectrometry*. ChemMedChem, 2006. **1**(7): p. 702-705.
4. Jones, R.W., R.B. Cody, and J.F. McClelland, *Differentiating Writing Inks Using Direct Analysis in Real Time Mass Spectrometry*. Journal of Forensic Sciences, 2006. **51**(4): p. 915-918.
5. Morlock, G. and W. Schwack, *Determination of isopropylthioxanthone (ITX) in milk, yoghurt and fat by HPTLC-FLD, HPTLC-ESI/MS and HPTLC-DART/MS*. Analytical and Bioanalytical Chemistry, 2006. **385**(3): p. 586-595.
6. Konuma, K., *Guide to Direct Analysis in Real Time (DART) Ionization Source*. Farumarushi, 2007. **43**(9): p. 903-905.
7. Kpegba, K., et al., *Analysis of Self-Assembled Monolayers on Gold Surfaces Using Direct Analysis in Real Time Mass Spectrometry*. Analytical Chemistry, 2007. **79**(14): p. 5479-5483.
8. Kusai, A., *Fundamental and application of the direct analysis in real time mass spectrometry*. Bunseki, 2007. **3**: p. 124-127.
9. Laramée, J.A. and R.B. Cody, *Chemi-ionization and Direct Analysis in Real Time (DART™) Mass Spectrometry*, in *The Encyclopedia of Mass Spectrometry Volume 6: Ionization Methods* M.L. Gross and R.M. Caprioli, Editors. 2007, Elsevier: Amsterdam. p. 377-387.
10. Laramée, J.A., et al., *Forensic Application of DART (Direct Analysis in Real Time) Mass Spectrometry*, in *Forensic Analysis on the Cutting Edge: New Methods for Trace Evidence Analysis*, R.D. Blackledge, Editor. 2007, Wiley-Interscience Hoboken, NJ.
11. Morlock, G. and Y. Ueda, *New coupling of planar chromatography with direct analysis in real time mass spectrometry* Journal of Chromatography A, 2007. **1143**(1-2): p. 243-251

12. Morlock, G. and Y. Ueda, *Coupling Planar Chromatography with Time-of-Flight Mass Spectrometry Using an Open-Air Ion Source*. LCGC: The Peak, 2007: p. 7-13.
13. Pierce, C.Y., et al., *Ambient generation of fatty acid methyl ester ions from bacterial whole cells by direct analysis in real time (DART) mass spectrometry*. Chemical Communications, 2007(8): p. 807 - 809.
14. Ropero-Miller, J.D., et al., *Comparison of the Novel Direct Analysis in Real Time Time-of-Flight Mass Spectrometry (AccuTOF-DART™) and Signature Analysis for the Identification of Constituents of Refined Illicit Cocaine*. Microgram Journal, 2007. 5(1-4): p. 5.
15. Saitoh, K., *Direct analysis for fragrance ingredients using DART-TOFMS*. Aroma Research, 2007. 8(4): p. 366-369.
16. Steiner, R.R., *AccuTOF-DART Mass Spectrometry*. 2007, Virginia Department of Forensic Science Controlled Substances Procedures Manual, Section 31.
17. Stout, P.R., et al., *Evaluation of urine samples utilizing direct analysis real time-of-flight mass spectrometry (AccuTOF-DART™) for postmortem toxicology screening.*, in *Society of Forensic Toxicologists (SOFT) Annual Meeting*. 2007: Raleigh-Durham, NC.
18. Vail, T.M., P.R. Jones, and O.D. Sparkman, *Rapid and Unambiguous Identification of Melamine in Contaminated Pet Food Based on Mass Spectrometry with Four Degrees of Confirmation*. Journal of Analytical Toxicology, 2007. 31(6): p. 304-312.
19. Alpmann, A. and G. Morlock, *Rapid and sensitive determination of acrylamide in drinking water by planar chromatography and fluorescence detection after derivatization with dansulfinic acid*. Journal of Separation Science, 2008. 31(1): p. 71-77.
20. Ayers, S., et al., *Pharmacokinetic analysis of anti-allergy and anti-inflammation bioactives in a nettle (Urtica dioica) extract*. Online Journal of Pharmacology and Pharmacokinetics, 2008. 5: p. 6-21.
21. Banerjee, S., et al., *Expression of tropane alkaloids in the hairy root culture of **Atropa acuminata** substantiated by DART mass spectrometric technique*. Biomedical Chromatography, 2008. 22(8): p. 830-834.
22. Banerjee, S., et al., *Analysis of cell cultures of **Taxus wallichiana** using direct analysis in real-time mass spectrometric technique*. Biomedical Chromatography, 2008. 22(3): p. 250-253.

23. Cajka, T., J. Hajslova, and K. Mastovska, *Mass spectrometry and hyphenated instruments in food analysis*, in *Handbook of Food Analysis Instruments*, S. Ötleş, Editor. 2008, CRC Press, Taylor & Francis Group. p. 197-228.
24. Cajka, T., et al., *GC–TOF-MS and DART–TOF-MS: Challenges in the Analysis of Soft Drinks*. LC/GC Europe, 2008. **21**(5): p. 250-256.
25. Coates, C.M., et al., *Flammable Solvent Detection Directly from Common Household Materials Yields Differential Results: An Application of Direct Analysis in Real-Time Mass Spectrometry*. *Journal of Forensic Identification*, 2008. **58**(6): p. 624 -631.
26. Fernandez, F.M., M.D. Green, and P.N. Newton, *Prevalence and Detection of Counterfeit Pharmaceuticals: A Mini Review*. *Ind. Eng. Chem. Res.*, 2008. **47**(3): p. 585-590.
27. Grange, A.H., *An Inexpensive Autosampler to Maximize Throughput for an Ion Source that Samples Surfaces in Open Air*. *Environmental Forensics*, 2008. **9**(2): p. 127-136.
28. Grange, A.H., *An Integrated Wipe Sample Transport/Autosampler to Maximize Throughput for a Direct Analysis in Real Time (DART)/Orthogonal Acceleration, Time-of-Flight Mass Spectrometer (oa-TOFMS)*. *Environmental Forensics*, 2008. **9**(2): p. 137 - 143.
29. Grange, A.H., *An Autosampler and Field Sample Carrier for Maximizing Throughput Using an Open-Air Source for MS*. *American Laboratory*, 2008. **40**(16): p. 11-13.
30. Grange, A.H. and G.W. Sovocool, *Automated determination of precursor ion, product ion, and neutral loss compositions and deconvolution of composite mass spectra using ion correlation based on exact masses and relative isotopic abundances*. *Rapid Communications in Mass Spectrometry*, 2008. **22**(15): p. 2375-2390.
31. Hajslova, J., et al. (2008) *Analysis of Deoxynivalenol in Beer*. *LCGC: Chromatography Online*.
32. Hajslova, J., et al., *DART–TOFMS: A Challenging Approach in Rapid Monitoring of Brominated Flame Retardants in Environmental Matrices*. *Organohalogen Compounds*, 2008. **70**: p. 922-925.
33. Harris, G.A., L. Nyadong, and F.M. Fernandez, *Recent developments in ambient ionization techniques for analytical mass spectrometry*. *Analyst*, 2008. **133**: p. 1297-1301.

34. Laramée, J.A., et al., *Detection of Chemical Warfare Agents on Surfaces Relevant to Homeland Security by Direct Analysis in Real-Time Spectrometry*. American Laboratory, 2008. **40**: p. 16-20.
35. Madhusudanan, K.P., et al., *Analysis of hairy root culture of Rauvolfia serpentina using direct analysis in real time mass spectrometric technique*. Biomedical Chromatography, 2008. **22**(6): p. 596-600.
36. Newton, P.N., et al., *A Collaborative Epidemiological Investigation into the Criminal Fake Artesunate Trade in South East Asia*. PLoS Medicine, 2008. **5**(2): p. e32.
37. Newton, P.N., et al., *Characterization of "Yaa Chud" Medicine on the Thailand–Myanmar Border: Selecting for Drug-resistant Malaria and Threatening Public Health*. American Journal of Tropical Medicine and Hygiene, 2008. **79**(5): p. 662-669.
38. Roschek Jr., B. and R.S. Alberte, *Pharmacokinetics of Cyanidin and Anti-Influenza Phytonutrients in an Elder Berry Extract Determined by LC-MS and DART TOF-MS*. Online Journal of Pharmacology and Pharmacokinetics, 2008. **4**: p. 1-17.
39. Schurek, J., et al., *Control of Strobilurin Fungicides in Wheat Using Direct Analysis in Real Time Accurate Time-of-Flight and Desorption Electrospray Ionization Linear Ion Trap Mass Spectrometry*. Analytical Chemistry, 2008. **80**(24): p. 9567-9575.
40. Smith, N.J., M.A. Domin, and L.T. Scott, *HRMS Directly From TLC Slides. A Powerful Tool for Rapid Analysis of Organic Mixtures*. Organic Letters, 2008. **10**(16): p. 3493-3496.
41. Vaclavik, L., et al., *Direct analysis in real time–time-of-flight mass spectrometry: Analysis of pesticide residues and environmental contaminants*. Chemicke Listy, 2008. **102**: p. s324–s327.
42. Yang, X., X. Xu, and H.-F. Ji, *Solvent Effect on the Self-Assembled Structure of an Amphiphilic Perylene Diimide Derivative*. Journal of Physical Chemistry B, 2008. **112**: p. 7196–7202.
43. Yew, J.Y., R.B. Cody, and E.A. Kravitz, *Cuticular hydrocarbon analysis of an awake behaving fly using direct analysis in real-time time-of-flight mass spectrometry*. Proceedings of the National Academy of Sciences, 2008. **105**(20): p. 7135-7140.
44. Bennett, M.J. and R.R. Steiner, *Detection of GHB in various drink matrices via AccuTOF-DART*. Journal of Forensic Science, 2009. **54**(2): p. 370-375.

45. Cody, R.B., *The Observation of Molecular Ions and Analysis of Nonpolar Compounds with the Direct Analysis in Real Time Ion Source*. Analytical Chemistry, 2009. **81**(3): p. 1101-1107.
46. Curtis, M.E., et al., *Determination of the Presence or Absence of Sulfur Materials in Drywall Using Direct Analysis in Real Time in Conjunction with an Accurate-Mass Time-of-Flight Mass Spectrometer*. Journal of the American Society for Mass Spectrometry, 2009. **20**(11): p. 2082-2086.
47. Dane, A.J. and R.B. Cody *Using Solid Phase Microextraction with AccuTOF-DART™ for Fragrance Analysis*, in *LC/GC: The Application Notebook*. 2009.
48. Grange, A.H., *Rapid Semi-Quantitative Surface Mapping of Airborne-Dispersed Chemicals Using Mass Spectrometry*. Environmental Forensics, 2009. **10**(3): p. 183-195.
49. Harris, G.A. and F.M. Fernandez, *Simulations and Experimental Investigation of Atmospheric Transport in an Ambient Metastable-Induced Chemical Ionization Source*. Analytical Chemistry, 2009. **81**(1): p. 322-329.
50. Jagerdeo, E. and M. Abdel-Rehim, *Screening of Cocaine and Its Metabolites in Human Urine Samples by Direct Analysis in Real-Time Source Coupled to Time-of-Flight Mass Spectrometry After Online Preconcentration Utilizing Microextraction by Packed Sorbent*. Journal of the American Society for Mass Spectrometry, 2009. **20**(5): p. 891-899.
51. Jee, E.H., et al., *Detection of characterising compounds on TLC by DART-MS*. Planta Med, 2009. **75**(09): p. PG38.
52. Kawamura, M., R. Kikura-Hanajiri, and Y. Goda, *Simple and rapid screening for psychotropic natural products using Direct Analysis in Real Time (DART)-TOFMS*. Yakugaku Zasshi : Journal of the Pharmaceutical Society of Japan, 2009. **129**(6): p. 719-25.
53. Kim, H.J. and Y.P. Jang, *Direct analysis of curcumin in turmeric by DART-MS*. Phytochemical Analysis, 2009. **published online**.
54. Konuma, K., *Elementary Guide to Ionization Methods for Mass Spectrometry—Introduction of Direct Analysis in Real Time Mass Spectrometry*. Bunseki, 2009. **9**: p. 464-467.
55. Laramée, J.A., et al., *Detection of Peroxide and Tetrazine Explosives on Surfaces by Direct Analysis in Real Time Mass Spectrometry*. American Laboratory Online Edition, 2009. **2**(2): p. 1-5.
56. Laramée, J.A., et al., *Alcohols Can Now Be Analyzed by a Direct Analysis in Real-Time Method: Applications for Chemical Warfare Agent Synthesis*. American Laboratory, 2009. **41**(4): p. 24-27.

57. Laramée, J.A., et al., *An Improved Protocol for the Analysis of Alcohols by Direct Analysis in Real Time Mass Spectrometry*. American Laboratory, 2009. **41**(7): p. 25-27.
58. Maleknia, S.D., T.L. Bell, and M.A. Adam, *Eucalypt smoke and wildfires: Temperature dependent emissions of biogenic volatile organic compounds*. International Journal of Mass Spectrometry, 2009. **279**(2-3): p. 126-133.
59. Maleknia, S.D., et al., *Temperature-dependent release of volatile organic compounds of eucalypts by direct analysis in real time (DART) mass spectrometry*. Rapid Communications in Mass Spectrometry, 2009. **23**(15): p. 2241-2246.
60. Mayoral, J.G., et al., *NADP<sup>+</sup>-dependent farnesol dehydrogenase, a corpora allata enzyme involved in juvenile hormone synthesis*. Proceedings of the National Academy of Sciences, 2009. **106**(50): p. 21091-21096.
61. Nilles, J.M., T.R. Connell, and H.D. Durst, *Quantitation of Chemical Warfare Agents Using the Direct Analysis in Real Time (DART) Technique*. Analytical Chemistry, 2009. **81**(16): p. 6744-6749.
62. Nyadong, L., et al., *Combining Two-Dimensional Diffusion-Ordered Nuclear Magnetic Resonance Spectroscopy, Imaging Desorption Electrospray Ionization Mass Spectrometry, and Direct Analysis in Real-Time Mass Spectrometry for the Integral Investigation of Counterfeit Pharmaceuticals*. Analytical Chemistry, 2009. **81**(12): p. 4803-4812.
63. Roschek Jr., B., et al., *Pro-Inflammatory Enzymes, Cyclooxygenase 1, Cyclooxygenase 2, and 5-Lipoxygenase, Inhibited by Stabilized Rice Bran Extracts*. Journal of Medicinal Food, 2009. **12**(3): p. 615–623.
64. Roschek Jr., B., et al., *Elderberry flavonoids bind to and prevent H1N1 infection in vitro* Phytochemistry, 2009. **70**(10): p. 1255-1262.
65. Saka, K., et al., *Identification of active ingredients in dietary supplements using non-destructive mass spectrometry and liquid chromatography–mass spectrometry*. Forensic Science International, 2009. **191**(1-3): p. e5-e10.
66. Song, L., et al., *Ionization Mechanism of Negative Ion-Direct Analysis in Real Time: A Comparative Study with Negative Ion-Atmospheric Pressure Photoionization*. Journal of the American Society for Mass Spectrometry, 2009. **20**(1): p. 42-50.
67. Song, L., et al., *Ionization Mechanism of Positive-Ion Direct Analysis in Real Time: A Transient Microenvironment Concept*. Analytical Chemistry, 2009. **81**(24): p. 10080-10088.

68. Steiner, R.R. and R.L. Larson, *Validation of the Direct Analysis in Real Time (DART) Source for Use in Forensic Drug Screening*. Journal of Forensic Science, 2009. **54**(3): p. 617-622.
69. Uchiyama, N., et al., *Identification of a cannabimimetic indole as a designer drug in a herbal product*. Forensic Toxicology, 2009. **27**(2): p. 61-66.
70. Vaclavik, L., et al., *Ambient mass spectrometry employing direct analysis in real time (DART) ion source for olive oil quality and authenticity assessment*. Analytica Chimica Acta, 2009. **645**(1-2): p. 56-63.
71. Wang, C., Y.-I. Zhou, and D.C. Baker, *Synthesis of a mannose-capped disaccharide with a thiol terminus* ARKIVOC, 2009. **xiv**: p. 171-180.
72. Bajpai, V., et al., *Profiling of Piper betle Linn. cultivars by direct analysis in real time mass spectrometric technique*. Biomedical Chromatography, 2010. **24**(12): p. 1283-1286.
73. Bevilacqua, V.L.H., et al., *Ricin Activity Assay by Direct Analysis in Real Time Mass Spectrometry Detection of Adenine Release*. Analytical Chemistry, 2010. **82**(3): p. 798-800.
74. Block, E., *Garlic and Other Alliums. The Lore and the Science*. 2010, Cambridge, UK: RSC Publishing.
75. Block, E., et al., *Allium chemistry: Use of new instrumental techniques to "see" reactive organosulfur species formed upon crushing garlic and onion*. Pure and Applied Chemistry, 2010. **82**(3): p. 535-539.
76. Block, E., et al., *Applications of Direct Analysis in Real Time Mass Spectrometry (DART-MS) in Allium Chemistry. 2-Propenesulfenic and 2-Propenesulfinic Acids, Diallyl Trisulfane S-Oxide, and Other Reactive Sulfur Compounds from Crushed Garlic and Other Alliums*. Journal of Agricultural and Food Chemistry, 2010. **58**(8): p. 4617-4625.
77. Cajka, T., et al., *Recognition of beer brand based on multivariate analysis of volatile fingerprint*. Journal of Chromatography A, 2010. **1217**(25): p. 4195-4203.
78. Cajka, T., et al., *Ambient mass spectrometry employing a DART ion source for metabolomic fingerprinting/profiling: A powerful tool for beer origin recognition*. Metabolomics, 2010. **published online: 9 December 2010**.
79. Chernetsova, E.S., et al., *An ultra superfast identification of low-molecular components of pharmaceuticals by DART mass spectrometry*. Journal of Analytical Chemistry, 2010. **65**(14): p. 1537-1539.

80. Chernetsova, E.S., et al., *Capabilities of direct analysis in real time mass spectrometry and gas chromatography-mass spectrometry in the mint oil test*. Mendeleev Communications, 2010. **20**(5): p. 299-300.
81. Chernetsova, E.S., et al., *The use of DART mass spectrometry for express confirmation of empirical formulas of heterocyclic compounds*. Russian Chemical Bulletin, 2010. **59**(10): p. 2014-2015.
82. Cody, R.B. and A.J. Dane, *Direct Analysis in Real Time Ion Source*, in *Encyclopedia of Analytical Chemistry*, R.A. Meyers, Editor. 2010, John Wiley & Sons, Ltd.: Published online: December 15.
83. Curtis, M., et al., *Direct Analysis in Real Time (DART) Mass Spectrometry of Nucleotides and Nucleosides: Elucidation of a Novel Fragment  $[C_5H_5O]^+$  and Its In-Source Adducts*. Journal of the American Society for Mass Spectrometry, 2010. **21**(8): p. 1371-1381.
84. Dane, A.J. and R.B. Cody, *Selective ionization of melamine in powdered milk by using argon direct analysis in real time (DART) mass spectrometry*. Analyst, 2010. **135**(4): p. 696-699.
85. Domin, M.A., et al., *Routine analysis and characterization of highly insoluble polycyclic aromatic compounds by direct analysis in real time mass spectrometry (DART)*. Analyst, 2010. **135**(4): p. 700-704.
86. Eberherr, W., et al., *Investigations on the Coupling of High-Performance Liquid Chromatography to Direct Analysis in Real Time Mass Spectrometry*. Analytical Chemistry, 2010. **82**(13): p. 5792-5796.
87. Harris, G.A., et al., *Comparison of the Internal Energy Deposition of Direct Analysis in Real Time and Electrospray Ionization Time-of-Flight Mass Spectrometry*. Journal of the American Society for Mass Spectrometry, 2010. **21**(5): p. 855-863.
88. Haunschmidt, M., et al., *Determination of organic UV filters in water by stir bar sorptive extraction and direct analysis in real-time mass spectrometry*. Analytical and Bioanalytical Chemistry, 2010. **397**(1): p. 269-275.
89. Haunschmidt, M., et al., *Rapid identification of stabilisers in polypropylene using time-of-flight mass spectrometry and DART as ion source*. Analyst, 2010. **135**(1): p. 80-85.
90. Jones, R.W., T. Reinot, and J.F. McClelland, *Molecular Analysis of Primary Vapor and Char Products during Stepwise Pyrolysis of Poplar Biomass*. Energy & Fuels, 2010. **24**(9): p. 5199-5209.
91. Kim, H.J., et al., *Identification of marker compounds in herbal drugs on TLC with DART-MS*. Archives of Pharmacal Research, 2010. **33**(9): p. 1355-1359.

92. Kim, H.J., et al., *Quantitative analysis of major dibenzocyclooctane lignans in schisandrae fructus by online TLC-DART-MS*. *Phytochemical Analysis*, 2010. **22**(3): p. 258-262.
93. Kubec, R., et al., *Applications of Direct Analysis in Real Time Mass Spectrometry (DART-MS) in Allium Chemistry. (Z)-Butanethial S-Oxide and 1-Butenyl Thiosulfinates and Their S-(E)-1-Butenylcysteine S-Oxide Precursor from Allium sicutum*. *Journal of Agricultural and Food Chemistry*, 2010. **58**(2): p. 1121-1128.
94. Navare, A.T., et al., *Rapid direct analysis in real time (DART) mass spectrometric detection of juvenile hormone III and its terpene precursors*. *Analytical and Bioanalytical Chemistry*, 2010. **398**(7): p. 3005-3013.
95. Nilles, J.M., T.R. Connell, and H.D. Durst, *Thermal separation to facilitate Direct Analysis in Real Time (DART) of mixtures*. *Analyst*, 2010. **135**(5): p. 883-886.
96. Nilles, J.M., et al., *Explosives Detection Using Direct Analysis in Real Time (DART) Mass Spectrometry*. *Propellants, Explosives, Pyrotechnics*, 2010. **35**(5): p. 446-451.
97. Pérez, J.J., et al., *Transmission-mode direct analysis in real time and desorption electrospray ionization mass spectrometry of insecticide-treated bednets for malaria control*. *Analyst*, 2010. **135**(4): p. 712-719.
98. Ra, J., et al., *Bambusae Caulis in Taeniam extract reduces ovalbumin-induced airway inflammation and T helper 2 responses in mice*. *Journal of Ethnopharmacology*, 2010. **128**(1): p. 241-247.
99. Steiner, R.R., *A Rapid Technique for the Confirmation of Iodine and Red Phosphorus Using Direct Analysis in Real Time and Accurate Mass Spectrometry*. *Microgram Journal*, 2010. **7**(1): p. 3-6.
100. Vaclavik, L., et al., *Rapid determination of melamine and cyanuric acid in milk powder using direct analysis in real time-time-of-flight mass spectrometry*. *Journal of Chromatography A*, 2010. **1217**(25): p. 4204-4211.
101. Zhou, M., et al., *Rapid Mass Spectrometric Metabolic Profiling of Blood Sera Detects Ovarian Cancer with High Accuracy*. *Cancer Epidemiology, Biomarkers & Prevention*, 2010.
102. Zhou, M., J.F. McDonald, and F.M. Fernández, *Optimization of a Direct Analysis in Real Time/Time-of-Flight Mass Spectrometry Method for Rapid Serum Metabolomic Fingerprinting*. *Journal of the American Society for Mass Spectrometry*, 2010. **21**(1): p. 68-75.
103. Adams, J., *Analysis of printing and writing papers by using direct analysis in real time mass spectrometry*. *International Journal of Mass Spectrometry*, 2011. **301**(1-3): p. 109-126.

104. Beißmann, S., et al., *High-performance liquid chromatography coupled to direct analysis in real time mass spectrometry: Investigations on gradient elution and influence of complex matrices on signal intensities*. Journal of Chromatography A, 2011. **1218**(31): p. 5180-5186.
105. Biniecka, M. and S. Caroli, *Analytical methods for the quantification of volatile aromatic compounds*. TrAC Trends in Analytical Chemistry, 2011. **30**(11): p. 1756-1770.
106. Block, E., *Challenges and Artifact Concerns in Analysis of Volatile Sulfur Compounds*, in *Volatile Sulfur Compounds in Food, ACS Symposium Series*. 2011, American Chemical Society. p. 35-63.
107. Block, E., A.J. Dane, and R.B. Cody, *Crushing Garlic and Slicing Onions: Detection of Sulfenic Acids and Other Reactive Organosulfur Intermediates from Garlic and Other Alliums using Direct Analysis in Real-Time Mass Spectrometry (DART-MS)*. Phosphorus, Sulfur, and Silicon and the Related Elements, 2011. **186**(5): p. 1085-1093.
108. Cajka, T., et al., *Ambient mass spectrometry employing a DART ion source for metabolomic fingerprinting/profiling: a powerful tool for beer origin recognition*. Metabolomics, 2011. **7**(4): p. 500-508.
109. Cajka, T., et al., *Direct analysis of dithiocarbamate fungicides in fruit by ambient mass spectrometry*. Food additives & contaminants. Part A, Chemistry, analysis, control, exposure & risk assessment, 2011. **28**(10): p. 1372-82.
110. Chernetsova, E.S., et al., *New approach to detecting counterfeit drugs in tablets by DART mass spectrometry*. Pharmaceutical Chemistry Journal, 2011. **45**(5): p. 306-308.
111. Chernetsova, E.S., et al., *Clarification of the composition of  $[M+18]^+$  ions in DART mass spectra of polyethylene glycol using high-resolution mass spectrometry*. Journal of Analytical Chemistry, 2011. **66**(13): p. 1348-1351.
112. Cho, D.S., et al., *Evaluation of direct analysis in real time mass spectrometry for onsite monitoring of batch slurry reactions*. Rapid Communications in Mass Spectrometry, 2011. **25**(23): p. 3575-3580.
113. Deroo, C.S. and R.A. Armitage, *Direct Identification of Dyes in Textiles by Direct Analysis in Real Time-Time of Flight Mass Spectrometry*. Analytical Chemistry, 2011. **83**(18): p. 6924-6928.
114. Grange, A.H. and G.W. Sovocool, *Detection of illicit drugs on surfaces using direct analysis in real time (DART) time-of-flight mass spectrometry*. Rapid Communications in Mass Spectrometry, 2011. **25**(9): p. 1271-1281.

115. Hajslova, J., T. Cajka, and L. Vaclavik, *Challenging applications offered by direct analysis in real time (DART) in food-quality and safety analysis*. TrAC Trends in Analytical Chemistry, 2011. **30**(2): p. 204-218.
116. Haunschmidt, M., et al., *Identification and semi-quantitative analysis of parabens and UV filters in cosmetic products by direct-analysis-in-real-time mass spectrometry and gas chromatography with mass spectrometric detection*. Analytical Methods, 2011. **3**: p. 99-104.
117. Howlett, S.E. and R.R. Steiner, *Validation of Thin Layer Chromatography with AccuTOF-DART™ Detection for Forensic Drug Analysis*. Journal of Forensic Sciences, 2011. **56**(5): p. 1261-1267.
118. Kalachova, K., et al., *Simplified and rapid determination of polychlorinated biphenyls, polybrominated diphenyl ethers, and polycyclic aromatic hydrocarbons in fish and shrimps integrated into a single method*. Analytica Chimica Acta, 2011. **707**(1-2): p. 84-91.
119. Kawamura, M., R. Kikura-Hanajiri, and Y. Goda, *Simple and Rapid Screening for Methamphetamine and 3,4-Methylene-dioxymethamphetamine (MDMA) and Their Metabolites in Urine Using Direct Analysis in Real Time (DART)-TOFMS*. Yakugaku Zasshi : Journal of the Pharmaceutical Society of Japan, 2011. **131**(5): p. 827-833.
120. Kim, H.J., W.S. Baek, and Y.P. Jang, *Identification of ambiguous cubeb fruit by DART-MS-based fingerprinting combined with principal component analysis*. Food Chemistry, 2011. **129**(3): p. 1305-1310.
121. Kim, H.J., et al., *Quantitative analysis of major dibenzocyclooctane lignans in schisandrae fructus by online TLC-DART-MS*. Phytochemical Analysis, 2011. **22**(3): p. 258-262.
122. Kim, S.W., et al., *A rapid, simple method for the genetic discrimination of intact Arabidopsis thaliana mutant seeds using metabolic profiling by direct analysis in real-time mass spectrometry*. Plant Methods, 2011. **7**(1): p. 1-10.
123. Kpegba, K., et al., *Epiatzelechin from the Root Bark of Cassia sieberiana: Detection by DART Mass Spectrometry, Spectroscopic Characterization, and Antioxidant Properties*. Journal of Natural Products, 2011. **74**(3): p. 455-459.
124. Kubec, R., et al., *Precursors and Formation of Pyrithione and Other Pyridyl-Containing Sulfur Compounds in Drumstick Onion, Allium stipitatum*. Journal of Agricultural and Food Chemistry, 2011. **59**(10): p. 5763-5770.
125. Kucerová, P., et al., *Allium discoloration: the precursor and formation of the red pigment in giant onion (Allium giganteum Regel) and some other subgenus Melanocrommyum species*. Journal of Agricultural and Food Chemistry, 2011. **59**(5): p. 1821-1828.

126. Lee, J.H., et al., *Danshen extract does not alter pharmacokinetics of docetaxel and clopidogrel, reflecting its negligible potential in P-glycoprotein- and cytochrome P4503A-mediated herb–drug interactions*. International Journal of Pharmaceutics, 2011. **410**(1-2): p. 68-74.
127. Mess, A., et al., *Qualitative Analysis of Tackifier Resins in Pressure Sensitive Adhesives Using Direct Analysis in Real Time Time-of-Flight Mass Spectrometry*. Analytical Chemistry, 2011. **83**(19): p. 7323-7330.
128. Pfaff, A.M. and R.R. Steiner, *Development and validation of AccuTOF-DART™ as a screening method for analysis of bank security device and pepper spray components*. Forensic Science International, 2011. **206**(1-3): p. 62-70.
129. Samms, W.C., et al., *Analysis of Alprazolam by DART-TOF Mass Spectrometry in Counterfeit and Routine Drug Identification Cases*. Journal of Forensic Sciences, 2011. **56**(4): p. 993-998.
130. Sanchez, L.M., et al., *Versatile Method for the Detection of Covalently Bound Substrates on Solid Supports by DART Mass Spectrometry*. Organic Letters, 2011. **13**(15): p. 3770-3773.
131. Vaclavik, L., et al., *Authentication of Animal Fats Using Direct Analysis in Real Time (DART) Ionization-Mass Spectrometry and Chemometric Tools*. Journal of Agricultural and Food Chemistry, 2011. **59**(11): p. 5919-5926.
132. Wood, J.L. and R.R. Steiner, *Purification of pharmaceutical preparations using thin-layer chromatography to obtain mass spectra with Direct Analysis in Real Time and accurate mass spectrometry*. Drug Testing and Analysis, 2011. **3**(6): p. 345-351.
133. Cody, R.B., et al., *Rapid Classification of White Oak (Quercus alba) and Northern Red Oak (Quercus rubra) by Using Pyrolysis Direct Analysis in Real Time (DART™) and Time-of-Flight Mass Spectrometry*. Journal of Analytical and Applied Pyrolysis, 2012. **95**: p. 134-137.
134. Dunham, S.J.B., P.D. Hooker, and R.M. Hyde, *Identification, extraction and quantification of the synthetic cannabinoid JWH-018 from commercially available herbal marijuana alternatives*. Forensic Science International, 2012. **223**(1-3): p. 241-244.
135. Grabenauer, M., et al., *Analysis of Synthetic Cannabinoids Using High-Resolution Mass Spectrometry and Mass Defect Filtering: Implications for Nontargeted Screening of Designer Drugs*. Analytical Chemistry, 2012. **84**(13): p. 5574-5581.
136. He, X.N., et al., *Mass spectrometry of solid samples in open air using combined laser ionization and ambient metastable ionization*. Spectrochimica Acta Part B, 2012. **67**: p. 64-73.

137. Hu, B., et al., *Selective hydrogenation of CO<sub>2</sub> and CO to useful light olefins over octahedral molecular sieve manganese oxide supported iron catalysts*. Applied Catalysis B: Environmental, 2012. **132-133**: p. 54-61.
138. Lancaster, C. and E. Espinoza, *Analysis of select Dalbergia and trade timber using direct analysis in real time and time-of-flight mass spectrometry for CITES enforcement*. Rapid Communications in Mass Spectrometry, 2012. **26**(9): p. 1147-1156.
139. Lee, S.M., H.-J. Kim, and Y.P. Jang, *Chemometric Classification of Morphologically Similar Umbelliferae Medicinal Herbs by DART-TOF-MS Fingerprint*. Phytochemical Analysis, 2012. **23**(5): p. 508-512.
140. Li, Y., *Confined direct analysis in real time ion source and its applications in analysis of volatile organic compounds of Citrus limon (lemon) and Allium cepa (onion)*. Rapid Communications in Mass Spectrometry, 2012. **26**(10): p. 1194-1202.
141. Musah, R.A., et al., *Direct Analysis in Real Time Mass Spectrometry for Analysis of Sexual Assault Evidence*. Rapid Communications in Mass Spectrometry, 2012. **26**(9): p. 1039-1046.
142. Musah, R.A., et al., *Direct analysis in real time mass spectrometry with collision-induced dissociation for structural analysis of synthetic cannabinoids*. Rapid Communications in Mass Spectrometry, 2012. **26**(19): p. 2335-2342.
143. Musah, R.A., et al., *Rapid identification of synthetic cannabinoids in herbal samples via direct analysis in real time mass spectrometry*. Rapid Communications in Mass Spectrometry, 2012. **26**(9): p. 1109-1114.
144. Novotná, H., et al., *Metabolomic fingerprinting employing DART-TOFMS for authentication of tomatoes and peppers from organic and conventional farming*. Food Additives & Contaminants: Part A, 2012. **29**(9): p. 1335-1346.
145. Singh, S. and S.K. Verma, *Study of the distribution profile of piperidine alkaloids in various parts of Prosopis juliflora by the application of Direct Analysis in Real Time Mass Spectrometry (DART-MS)*. Natural Products and Bioprospecting, 2012. **2**(5): p. 206-209.
146. Singh, V., et al., *Direct Analysis in Real Time by Mass Spectrometric Technique for Determining the Variation in Metabolite Profiles of Cinnamomum tamala Nees and Eberm Genotypes*. The Scientific World Journal, 2012. **2012**: p. 6.
147. Uchiyama, N., et al., *URB-754: A new class of designer drug and 12 synthetic cannabinoids detected in illegal products*. Forensic Science International, 2012. **Publication Date (Web): 9 October 2012**.

148. Bentayeb, K., et al., *Non-visible print set-off of photoinitiators in food packaging: detection by ambient ionisation mass spectrometry*. Food Additives & Contaminants: Part A, 2013. **30**(4): p. 750-759.
149. Cajka, T., et al., *Evaluation of direct analysis in real time ionization-mass spectrometry (DART-MS) in fish metabolomics aimed to assess the response to dietary supplementation*. Talanta, 2013. **115**: p. 263-270.
150. Cajka, T., et al., *Application of direct analysis in real time ionization-mass spectrometry (DART-MS) in chicken meat metabolomics aiming at the retrospective control of feed fraud*. Metabolomics, 2013. **9**(3): p. 545-557.
151. Cody, R.B., *What Is the Opposite of Pandora's Box? Direct Analysis, Ambient Ionization, and a New Generation of Atmospheric Pressure Ion Sources*. Mass Spectrometry, 2013. **2**(Special Issue): p. S0007.
152. Cody, R.B. and A.J. Dane, *Soft Ionization of Saturated Hydrocarbons, Alcohols and Nonpolar Compounds by Negative-Ion Direct Analysis in Real-Time Mass Spectrometry*. Journal of the American Society for Mass Spectrometry, 2013. **24**(3): p. 329-334.
153. Djelal, H., et al., *The use of HPTLC and Direct Analysis in Real Time-Of-Flight Mass Spectrometry (DART-TOF-MS) for rapid analysis of degradation by oxidation and sonication of an azo dye*. Arabian Journal of Chemistry, 2013. **published online: 13 June 2013**.
154. Dwivedi, P., et al., *Electro-Thermal Vaporization Direct Analysis in Real Time-Mass Spectrometry for Water Contaminant Analysis during Space Missions*. Analytical Chemistry, 2013. **85**(20): p. 9898-9906.
155. Fukuda, E., et al., *Application of mixture analysis to crude materials from natural resources (IV)[1(a-c)]: identification of Glycyrrhiza species by direct analysis in real time mass spectrometry (II)*. Natural Product Communications, 2013. **8**(12): p. 1721-1724.
156. Grange, A.H., *Semi-quantitative analysis of contaminants in soils by direct analysis in real time (DART) mass spectrometry*. Rapid Communications in Mass Spectrometry, 2013. **27**(2): p. 305-318.
157. Hintersteiner, I., R. Hertsens, and C.W. Klampfl, *Direct Analysis in Real Time/Time-of-Flight Mass Spectrometry: Investigations on Parameters for the Coupling with Liquid-Phase Sample Introduction Techniques*. Journal of Liquid Chromatography & Related Technologies, 2013. **37**(13): p. 1862-1872.
158. Houlgrave, S., et al., *The Classification of Inkjet Inks Using AccuTOF™DART™ (Direct Analysis in Real Time) Mass Spectrometry—A Preliminary Study*. Journal of Forensic Sciences, 2013. **58**(3): p. 813-821.

159. Jones, R.W. and J.F. McClelland, *Analysis of writing inks on paper using direct analysis in real time mass spectrometry*. Forensic Science International, 2013. **231**(1-3): p. 73-81.
160. Lesiak, A.D., et al., *Direct Analysis in Real Time Mass Spectrometry (DART-MS) of "Bath Salt" Cathinone Drug Mixtures*. Analyst, 2013. **138**(12): p. 3424-3432.
161. Li, Y., *Applications of a confined DART (direct analysis in real time) ion source for online in vivo analysis of human breath*. Analytical Methods, 2013. **5**(24): p. 6933-6940.
162. Lim, A.Y., et al., *Detection of drugs in latent fingerprints by mass spectrometric methods*. Analytical Methods, 2013. **5**(17): p. 4378-4385.
163. Mess, A., et al., *A novel sampling method for identification of endogenous skin surface compounds by use of DART-MS and MALDI-MS*. Talanta, 2013. **103**: p. 398-402.
164. Park, H.M., et al., *Direct Analysis in Real Time Mass Spectrometry (DART-MS) Analysis of Skin Metabolome Changes in the Ultraviolet B-Induced Mice*. Biomolecules and Therapeutics, 2013. **21**(6): p. 470-475.
165. Rajchl, A., et al., *Rapid determination of 5-hydroxymethylfurfural by DART ionization with time-of-flight mass spectrometry*. Analytical and Bioanalytical Chemistry, 2013. **405**(14): p. 4737-4745.
166. Sisco, E., J. Dake, and C. Bridge, *Screening for trace explosives by AccuTOF™-DART®: An in-depth validation study*. Forensic Science International, 2013. **232**(1-3): p. 160-168.
167. Swider, J.R., *Optimizing Accu Time-of-Flight/Direct Analysis in Real Time for Explosive Residue Analysis*. Journal of Forensic Sciences, 2013. **58**(6): p. 1601-1606.
168. Al-Balaa, D., et al., *DART mass spectrometry for rapid screening and quantitative determination of cholesterol in egg pasta*. Journal of Mass Spectrometry, 2014. **49**(9): p. 911-917.
169. Easter, J.L. and R.R. Steiner, *Pharmaceutical identifier confirmation via DART-TOF*. Forensic Science International, 2014. **240**: p. 9-20.
170. Espinoza, E.O., et al., *Distinguishing Wild from Cultivated Agarwood (Aquilaria spp.) Using Direct Analysis in Real Time (DART™) and Time-of-Flight Mass Spectrometry*. Rapid Communications in Mass Spectrometry, 2014. **28**(3): p. 281-289.

171. Fukuda, E., et al., *Identification of the country of growth of Sophora flavescens using direct analysis in real time mass spectrometry (DART-MS)*. Natural Product Communications, 2014. **9**(11): p. 1591-4.
172. Jacobs, A.D. and R.R. Steiner, *Detection of the Duquenois–Levine chromophore in a marijuana sample*. Forensic Science International, 2014. **239**: p. 1-5.
173. Kiguchi, O., et al., *Thin-layer chromatography/direct analysis in real time time-of-flight mass spectrometry and isotope dilution to analyze organophosphorus insecticides in fatty foods*. Journal of Chromatography A, 2014. **1370**: p. 246-254.
174. Kim, H., et al., *DART-TOF-MS based metabolomics study for the discrimination analysis of geographical origin of Angelica gigas roots collected from Korea and China*. Metabolomics, 2014: p. 1-7.
175. Kim, H.J., S.R. Park, and Y.P. Jang, *Extraction-free In situ Derivatisation of Timosaponin AIII Using Direct Analysis in Real Time TOF/MS*. Phytochemical Analysis, 2014. **25**(4): p. 373-377.
176. Lesiak, A.D., et al., *DART-MS for rapid, preliminary screening of urine for DMAA*. Drug Testing and Analysis, 2014. **6**(7-8): p. 788-796.
177. Lesiak, A.D., et al., *Rapid detection by direct analysis in real time-mass spectrometry (DART-MS) of psychoactive plant drugs of abuse: The case of Mitragyna speciosa aka "Kratom"*. Forensic Science International, 2014. **242**: p. 210-218.
178. Lesiak, A.D., et al., *DART-MS as a Preliminary Screening Method for "Herbal Incense": Chemical Analysis of Synthetic Cannabinoids*. Journal of Forensic Science, 2014. **59**(2): p. 337-343.
179. Lesiak, A.D. and J.R.E. Shepard, *Recent advances in forensic drug analysis by DART-MS*. Bioanalysis, 2014. **6**(6): p. 819-842.
180. Morlock, G.E., P. Ristivojevic, and E.S. Chernetsova, *Combined multivariate data analysis of high-performance thin-layer chromatography fingerprints and direct analysis in real time mass spectra for profiling of natural products like propolis*. Journal of Chromatography A, 2014. **1328**: p. 104-112.
181. Musah, R.A., et al., *DART-MS in-source collision induced dissociation and high mass accuracy for new psychoactive substance determinations*. Forensic Science International, 2014. **244**: p. 42-49.
182. Paseiro-Cerrato, R., G.O. Noonan, and T.H. Begley, *Development of a rapid screening method to determine primary aromatic amines in kitchen utensils using direct analysis in real time mass spectrometry (DART-MS)*. Food Additives & Contaminants: Part A, 2014. **31**(3): p. 537-545.

183. Prchalová, J., et al., *Characterization of mustard seeds and paste by DART ionization with time-of-flight mass spectrometry*. Journal of Mass Spectrometry, 2014. **49**(9): p. 811-818.
184. Alusta, P., et al., *Instrumental improvements and sample preparations that enable reproducible, reliable acquisition of mass spectra from whole bacterial cells*. Rapid Communications in Mass Spectrometry, 2015. **29**(21): p. 1961-1968.
185. Armitage, R.A., C. Day, and K.A. Jakes, *Identification of anthraquinone dye colourants in red fibres from an Ohio Hopewell mound site by direct analysis in real time mass spectrometry*. STAR: Science & Technology of Archaeological Research, 2015. **1**(2): p. 1-10.
186. Bag, A. and R.R. Chattopadhyay, *Evaluation of Synergistic Antibacterial and Antioxidant Efficacy of Essential Oils of Spices and Herbs in Combination*. PLoS ONE, 2015. **10**(7): p. e0131321.
187. Cody, R.B., et al., *Identification of bacteria by fatty acid profiling with direct analysis in real time mass spectrometry*. Rapid Communications in Mass Spectrometry, 2015. **29**(21): p. 2007-2012.
188. Curtis, M., et al., *Schlieren visualization of fluid dynamics effects in direct analysis in real time mass spectrometry*. Rapid Communications in Mass Spectrometry, 2015. **29**(5): p. 431-439.
189. Forbes, T.P. and E. Sisco, *Trace detection and competitive ionization of erythritol tetranitrate in mixtures using direct analysis in real time mass spectrometry*. Analytical Methods, 2015. **7**(8): p. 3632-3636.
190. Kumar, S., et al., *Rapid fingerprinting of Rauwolfia species using direct analysis in real time mass spectrometry combined with principal component analysis for their discrimination*. Analytical Methods, 2015. **7**(14): p. 6021-6026.
191. Le Pogam, P., et al., *Analysis of the cyanolichen Lichina pygmaea metabolites using in situ DART-MS: from detection to thermochemistry of mycosporine serinol*. Journal of Mass Spectrometry, 2015. **50**(3): p. 454-462.
192. Lesiak, A.D., et al., *Plant Seed Species Identification from Chemical Fingerprints: A High-Throughput Application of Direct Analysis in Real Time Mass Spectrometry*. Analytical Chemistry, 2015. **87**(17): p. 8748-8757.
193. Lu, Y., et al., *Sensitivity and intensity enhancement in open air mass spectrometry assisted with a continuous wave infrared laser*. Journal of Analytical Atomic Spectrometry, 2015. **30**(7): p. 1663-1667.
194. Marino, M.A., et al., *Rapid Identification of Synthetic Cannabinoids in Herbal Incenses with DART-MS and NMR*. Journal of Forensic Sciences, 2015. **61**(S1): p. S82-S91.

195. McClure, P.J., G.D. Chavarria, and E. Espinoza, *Metabolic chemotypes of CITES protected Dalbergia timbers from Africa, Madagascar, and Asia*. Rapid Communications in Mass Spectrometry, 2015. **29**(9): p. 783-788.
196. Musah, R.A., et al., *A High Throughput Ambient Mass Spectrometric Approach to Species Identification and Classification from Chemical Fingerprint Signatures*. Scientific Reports, 2015. **5**: p. 11520\_1-16.
197. Poklis, J.L., et al., *Analysis of 25I-NBOMe, 25B-NBOMe, 25C-NBOMe and Other Dimethoxyphenyl-N-[(2-Methoxyphenyl) Methyl]Ethanamine Derivatives on Blotter Paper*. Journal of Analytical Toxicology, 2015. **39**(8): p. 617-623.
198. Poklis, J.L., et al., *Analysis of the First- and Second-Generation Raving Dragon Novelty Bath Salts Containing Methylone and Pentedrone*. Journal of Forensic Sciences, 2015. **60**(s1): p. S234-S240.
199. Rajchl, A., et al., *Evaluation of ice-tea quality by DART-TOF/MS*. Journal of Mass Spectrometry, 2015. **50**(11): p. 1214-1221.
200. Singh, A., et al., *Rapid screening and distribution of bioactive compounds in different parts of Berberis petiolaris using direct analysis in real time mass spectrometry*. Journal of Pharmaceutical Analysis, 2015. **5**(5): p. 332-335.
201. Sisco, E. and T.P. Forbes, *Rapid detection of sugar alcohol precursors and corresponding nitrate ester explosives using direct analysis in real time mass spectrometry*. Analyst, 2015. **140**(8): p. 2785-2796.
202. Song, Y., et al., *Rapidly differentiating grape seeds from different sources based on characteristic fingerprints using direct analysis in real time coupled with time-of-flight mass spectrometry combined with chemometrics*. Journal of Separation Science, 2015. **38**(17): p. 3084-3089.
203. Yew, J.Y. and H. Chung, *Insect pheromones: An overview of function, form, and discovery*. Progress in Lipid Research, 2015. **59**: p. 88-105.
204. Zhou, S., M.W. Forbes, and J.P.D. Abbatt, *Application of Direct Analysis in Real Time-Mass Spectrometry (DART-MS) to the Study of Gas-Surface Heterogeneous Reactions: Focus on Ozone and PAHs*. Analytical Chemistry, 2015. **87**(9): p. 4733-4740.
205. Cody, R.B. and A.J. Dane, *Alternative mass reference standards for direct analysis in real time mass spectrometry*. Rapid Communications in Mass Spectrometry, 2016. **30**(10): p. 1206-1212.
206. Cody, R.B. and A.J. Dane, *Dopant-assisted direct analysis in real time mass spectrometry with argon gas*. Rapid Communications in Mass Spectrometry, 2016. **30**(10): p. 1181-1189.

207. Lesiak, A.D., et al., *Direct analysis in real time high resolution mass spectrometry as a tool for rapid characterization of mind-altering plant materials and revelation of supplement adulteration - The case of Kanna*. Forensic Science International, 2016. **260**: p. 66-73.
208. Lesiak, A.D. and R.A. Musah, *More than just heat: ambient ionization mass spectrometry for determination of the species of origin of processed commercial products-application to psychoactive pepper supplements*. Analytical Methods, 2016. **8**(7): p. 1646-1658.
209. Maric, M. and C. Bridge, *Characterizing and Classifying Water-Based Lubricants using Direct Analysis in Real Time<sup>®</sup> - Time of Flight Mass Spectrometry*. Forensic Science International, 2016. **published online: 11 May 2016**.
210. Peace, M.R., et al., *Concentration of Nicotine and Glycols in 27 Electronic Cigarette Formulations*. Journal of Analytical Toxicology, 2016. **published online: May 10, 2016**.
211. Peace, M.R., et al., *Analysis of a Commercial Marijuana e-Cigarette Formulation*. Journal of Analytical Toxicology, 2016. **published online: April 7, 2016**.
212. Roberts, E.S., et al., *Analysis of carbohydrates in Fusarium verticillioides using size-exclusion HPLC - DRI and direct analysis in real time ionization - time-of-flight - mass spectrometry (DART-MS)*. Analytical Methods, 2016. **8**(3): p. 673-681.
213. Sisco, E. and J. Dake, *Detection of low molecular weight adulterants in beverages by direct analysis in real time mass spectrometry*. Analytical Methods, 2016. **8**(14): p. 2971-2978.
214. Sisco, E. and T.P. Forbes, *Direct analysis in real time mass spectrometry of potential by-products from homemade nitrate ester explosive synthesis*. Talanta, 2016. **150**: p. 177-183.

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Copyright © 2007 JEOL Ltd.