Solutions for Innovation

JMS-S3000 SpiralTOF

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
Fig. 1 Reduced topographic effect of matrix crystal.

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Table 1. Sample information and preparation conditions.

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
Fig. 4 MALDI mass spectra of PEG1500


Fig. 5 Averaged mass resolving power (10 spots) for  $m/z$  1097.6,  $m/z$  1537.9 and  $m/z$  1978.2.


Fig. 6 Mass error with external calibration method at  $m/z$  1097.6,  $m/z$  1537.9 and  $m/z$  1978.2 in eight spots mass spectra using four typical matrix compounds.

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
Fig. 7 SEM images of CHCA crystal with PEG1500: left: x500, right: x2,000.


Fig. 8 SEM images of DHB crystal with PEG1500: left: x500, right: x2,000



Fig. 9 SEM images of DIT crystal with PEG1500: left: x500, right: x2,000.




Fig. 10 SEM images of IAA crystal with PEG1500: left: x500, right: x2,000.

# Structural Analysis of Triacylglycerols by Using a MALDI-TOF/TOF System with Monoisotopic Precursor Selection

Ayumi Kubo,<sup>1</sup> Takaya Satoh,<sup>1</sup> Yoshiyuki Itoh,<sup>1</sup> Masahiro Hashimoto,<sup>1</sup> Jun Tamura,<sup>1</sup> Robert B. Cody<sup>2</sup>

<sup>1</sup>JEOL Ltd., Tokyo, Japan

<sup>2</sup>JEOL USA, Inc., Peabody, MA 01960, USA

**Abstract.** A new MALDI-TOF/TOF system with monoisotopic precursor selection was applied to the analysis of triacylglycerols in an olive oil sample. Monoisotopic precursor selection made it possible to obtain production mass spectra without interference from species that differed by a single double bond. Complete structure determination of all triacylglycerols, including structural isomers, was made possible by interpreting the charge-remote fragmentation resulting from high-energy collision-induced dissociation (CID) of the sodiated triacylglycerols.

**Key words:** MALDI, TOF/TOF, Charge-remote fragmentation, High-energy CID, Triacylglycerols, MS/MS, Lipids, Collision-induced dissociation

Received: 29 March 2012/Revised: 22 June 2012/Accepted: 22 June 2012/Published online: 18 December 2012

## Introduction

Triacylglycerols (TAGs or triglycerides) are comprised of three fatty acids esterified with glycerol. Because TAGs are the major components in animal fats and vegetable oils, the analysis of TAGs is biologically important and crucial for quality control of food products.

Recent mass spectrometric approaches to the analysis of TAGs have made use of atmospheric pressure chemical ionization (APCI) [1–5] or electrospray ionization (ESI) [6–10] and tandem mass spectrometry (MS/MS).

Using ESI and a triple quadrupole mass spectrometer, Hsu and Turk [6] reported that collisional activation of lithium adducts of TAGs can provide structural information about the acyl groups. Low-energy collision-induced dissociation (CID) of cationized TAGs does not provide information about the position of the double bonds. However, the CID fragments of unusual dilithiated species was shown to be dependent on double bond location. Byrdwell and Neff [7] reported a method based on dual parallel ESI and APCI combined with tandem mass spectrometry for the analysis of TAGs and their oxidation products. McAnoy et al. [8] used ESI with a linear ion trap to characterize TAG components within a complex mixture of neutral lipids from cell extracts.

High-energy CID is an especially attractive approach for TAG analysis because charge-remote fragmentation [11–22] provides a great deal of information about lipid structure. The complete structural characterization of TAGs was reported in 1998 by Cheng et al. using fast atom bombardment (FAB) and tandem magnetic sector mass spectrometry with high-energy CID fragmentation of the  $[M + Na]^+$  species [21]. All TAG structural features could be determined except stereochemistry.

However, large tandem magnetic sector mass spectrometers have fallen out of favor in recent years and high-energy CID appeared destined to become a “lost art” until the introduction of tandem time-of-flight (TOF/TOF) mass spectrometers by Cotter and Cornish in 1993 [22]. Recently, Pittenauer and Allmaier showed that TOF/TOF mass spectrometers have the potential to provide the same complete structural information as a tandem magnetic sector mass spectrometer [23]. The principal limitation of this method was found to be the poor MS-I selectivity (a 4 to 6 u window) of the TOF/TOF system, making it impractical to select precursor ions for TAGs with compositions that differ by two hydrogens. The authors concluded that a LC/MALDI-MS/MS approach might be required to make use of charge-remote fragmentations to characterize TAGs in complex mixtures.

We have developed a tandem time-of-flight mass spectrometer featuring high precursor ion selectivity that resolves the problem of poor MS-I selectivity [24]. The mass

Correspondence to: Robert Cody; e-mail: cody@jeol.com

spectrometer uses multi-turn and “perfect focusing” ion optics [25] to fit a very long (17-m) flight path into a compact space [26]. In TOF/TOF mode, an ion gate positioned at the 15 m point in the spiral ion flight path is used to isolate and guide the precursor ion into a gas-filled collision chamber. The long flight path provides ample time separation prior to precursor ion selection, resulting in unit precursor selectivity. The precursor ions undergo 20 kV collisions with a target gas and are subjected to a 9 kV post-acceleration into an offset parabolic reflectron with wide energy acceptance.

Monoisotopic precursor selection combined with high-energy CID is the key to using TOF/TOF for structural analysis of triacylglycerols in complex mixtures. This paper describes the method for structural analysis with this system and reports the complete structural analysis of TAGs, including isomers, in a commercial olive oil sample.

## Experimental

### Materials and Chemicals

A triacylglycerol standard (1-palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol), matrix (2,5-dihydroxybenzoic acid or DHB), and cationizing agent (sodium trifluoroacetate), were purchased from Sigma-Aldrich (St. Louis, MO). A trioleoyl-glycerol (triolein) standard was purchased from TCI. Tetrahydrofuran (THF) was purchased from Wako (Osaka, Japan) and olive oil was purchased from local stores. The triacylglycerol standard, including triolein, and the olive oil were dissolved in THF at respective concentrations of 100 pmol/uL and 10 ug/uL. A solution of sodium trifluoroacetate and DHB was dissolved in THF at respective concentrations of 1 ug/uL and 20 ug/uL, and added to the samples at a volume ratio of 1:1:2. The resulting mixture was loaded onto an MTP 96-hole hairline plate (JEOL Ltd., Akishima Japan) at a volume of 1 uL per spot.

### MALDI Mass Spectrometry

A JMS-S3000 Spiral TOF (JEOL Ltd., Akishima, Japan) equipped with the TOF/TOF option was used for all measurements. The laser was a Nd-YLF laser operated at

a wavelength of 349 nm. The laser intensity and the detector voltage were set to prevent triacylglycerol peaks from saturating. The extraction delay was optimized to 400 ns to provide a resolving power (FWHM) of approximately 50,000 for the TAG peaks in MS-I mode. For product-ion mass spectrum acquisition, helium collision gas was introduced to attenuate the precursor ion abundance to approximately 50 % of the initial value. The laser was operated at a repetition rate of 1000 Hz. Spectra were acquired at a rate of two spectra per s and 500 spectra were accumulated for each product-ion mass spectrum shown here.

## Results and Discussion

Figure 1 shows the structure of 1-palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol. The structure shows (18:1) oleic acid, (16:0) palmitic acid, and (18:2) linoleic acid substituents at position *sn-2* (the site that determines the stereochemistry) and positions *sn-1* and *sn-3*, respectively. In this article, we have labeled the fatty acid substituents at positions *sn-1* and *sn-3* as “*sn-1/sn-3*.” The substituents at *sn-1* and *sn-3* are indistinguishable by mass spectrometry because the steric structure of triacylglycerol cannot be identified by mass spectrometry. Each fragmentation path is assigned as shown in Figure 1, and is labeled alphabetically. Each letter represents the initial letter of the fatty acid, and the accompanying number represents the bonding position in each fatty acid. The labeling for TAGs such as “TAG(54:3)” follows the convention where the numeral on the left in parentheses represents the total number of acyl carbon chains and the numeral on the right represents the total number of unsaturated bonds at fatty acid moieties.

The major species observed for 1-palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol was the sodiated molecule  $[M + Na]^+$ . Figure 2 shows the product ion spectrum acquired by selecting the monoisotopic ion of this species. The resulting fragment ions are solely monoisotopic ions as well because a monoisotopic precursor ion was selected. Thus, each fragmentation path is observed as a single peak on the product-ion mass spectrum. Figure 2a shows the entire mass range of the product-ion spectrum. The  $Na^+$  peak detected at  $m/z$  23.0 confirms that the precursor ion is indeed  $[M + Na]^+$ . Peaks characteristic of fatty acid fragmentation are predicted as A-, B-,


Figure 1. Structure and charge-remote fragmentation of sodiated 1-palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol


Figure 2. Product-ion mass spectrum for sodiated 1-palmitoyl-2-oleoyl-3-linoleoyl-rac-glycerol, (a) entire mass range; (b)  $m/z$  650–890 magnified

C-, E-, G-, and J-type ions using the nomenclature defined in reference [21]. Figure 2a demonstrates that all of A-, B-, C-, E-, G-, and J-type ions predicted in reference [21] were observed for this example and “G+2” ions (mentioned in reference [21]) were observed. The structure of “G+2” ions and their fragmentation pathway are not clear, but “G+2” ions were also observed in the product ion spectrum of the of triolein standard (shown in Figure 5b) at a relatively lower intensity than that of the G-type ions. Figure 2b also shows that signals resulting from charge-remote fragmentation were detected in the high mass range above  $m/z$  650. When the fragment ion at each bonding position is defined as in Figure 1, the peaks can be assigned as shown in Figure 2b. The intensities of fragment ions corresponding to unsaturated bonding positions, such as L $\Delta$ 9, L $\Delta$ 12, and O $\Delta$ 9, are relatively weak or are not observed, resulting in a peak pattern that reflects the structure of 3 fatty acids.

In the analysis of triacylglycerols in the olive oil sample, particular attention was focused on the G- and J-type ions. These ions have the structure where two molecules of fatty acid are eliminated from the precursor ion [21]. These ions

help determine the numbers of carbon chains and unsaturated bonds in each fatty acid. In the G-type ion, fatty acids remain at *sn-1/sn-3*, while the J-type ion, where a fatty acid remains at *sn-2*, has one less CH<sub>2</sub> at the end. This makes it possible to estimate the bonding positions of three fatty acids because fatty acids having an odd acyl carbon number rarely exist in the natural world.

Figure 3 shows the mass spectrum of the olive oil. Sodiated triacylglycerols [M + Na]<sup>+</sup> were observed for this sample that included TAG (52:3) ( $m/z$  879.7), TAG (52:2) ( $m/z$  881.7), TAG (54:4) ( $m/z$  905.8), and TAG (54:3) ( $m/z$  907.8). The monoisotopic ions of these four TAGs were selected as the precursor ions, and their product-ion mass spectra were acquired. Figure 4 shows the spectra of ions at  $m/z$  905.8 acquired before and after the precursor ion selection. The figure demonstrates that only the ions at  $m/z$  905.8 were selected, completely eliminating ions at other mass values.

Figure 5 shows comparison between the product-ion mass spectra of TAG (54:3) [M + Na]<sup>+</sup> at  $m/z$  907.8 from olive oil and from triolein standard. Given that olive oil is


Figure 3. Mass spectrum of olive oil sample showing sodiated TAGs


Figure 4. Precursor-ion selection for ions at  $m/z$  905.8 (a) before selection, (b) after selection

rich in oleic acid, the ion at  $m/z$  907.8 is expected to contain three oleic acids (18:1). Both of the product-ion mass spectra show a  $J_2$ -type ion at  $m/z$  331.3, indicating that an oleic acid is bonded at the *sn*-2 position, and a G-type ion at  $m/z$  345.3, indicating that an oleic acid is bonded at the *sn*-1/*sn*-3 positions. The spectra show only one peak that is considered

an A-, B-, and C-type ion, suggesting that TAG (54:3) is trioleoylglycerol, which contains three oleic acid molecules. The signals in high mass region resulting from charge-remote fragmentation were identical between the sample and standard, and the spectral patterns were consistent with structure of oleic acid.


Figure 5. Comparison of product-ion mass spectra for the precursor at  $m/z$  907.8 with that of sodiated triolein standard, (a) the ions at  $m/z$  907.8 from olive oil, (b) triolein standard


Figure 6. Product-ion mass spectrum for the precursor at  $m/z$  905.8

Next, the ion at  $m/z$  905.8 was selected as the precursor ion. The  $m/z$  value of this ion suggests that it is a monoisotopic  $[M + Na]^+$  ion of TAG (54:4). Figure 6 shows the product-ion mass spectrum. It is expected that this triacylglycerol is also composed of two (18:1) oleic acids and one (18:2) linoleic acid, given that the major component of olive oil is oleic acid. The product-ion mass spectrum shows fragment ions assigned as G-type ions, at  $m/z$  343.4 and  $m/z$  345.3. If the ion at  $m/z$  343.4 is a G-type ion, the fatty acid molecule at  $sn-1/sn-3$  is linoleic acid, and if the ion at  $m/z$  345.4 is a G-type ion, the fatty acid molecules at  $sn-1/sn-3$  are oleic acid. The ion at  $m/z$  347.3 is assigned as a “G+2” ion because in the product-ion mass spectrum “G+2” ions were observed at lower intensity than G-type ions as discussed above, and the intensity of the ion at  $m/z$  347.3 is relatively lower than that of the ion at  $m/z$  345.3. This is consistent with the assignment of G+2 ions by Cheng et al. in reference [21]. Since the G-type ion suggests that both oleic acid and linoleic acid are bonded, the remaining fatty acid is (18:1) oleic acid. Next, the product-ion spectrum shows J-type ions: a J-type ion containing oleic acid ( $J_{2O}$ ) and a J-type ion containing linoleic acid ( $J_{2L}$ ) at  $m/z$  331.3 and  $m/z$  329.3, respectively. In the high-mass region, the signals resulting from charge-remote fragmentation were consistent with the structures of oleic acid and linoleic acid. This demonstrates that  $m/z$  905.8 is triacylglycerol that contains two molecules of oleic acid and one molecule of linoleic acid and is a mixture of the structural isomers 1,3-dioleoyl-2-linoleoyl-glycerol and 1,2-dioleoyl-3-linoleoyl-glycerol. Table 1 summarizes the structures of triacylglycerols determined for the olive oil samples from the product-ion mass spectra. The ions associated with the peak at  $m/z$  879.7 are a mixture of the structural isomers 1-palmitoyl-2-oleoyl-3-linoleoyl-glycerol and 1-palmitoyl-2-linoleoyl-3-oleoyl-glycerol.

Table 1. Summary of TAGs Found in the Olive Oil Sample

| $m/z$ | Acyl carbon number and number of double bond | Composition of each fatty acid    |
|-------|--|-----------------------------------|
| 879.7 | 52:3   | (16:0,18:1,18:2) (16:0,18:2,18:1) |
| 881.7 | 52:2   | (16:0,18:1,18:1)                  |
| 905.8 | 54:4   | (18:1,18:1,18:2) (18:1,18:2,18:1) |
| 907.8 | 54:3   | (18:1,18:1,18:1)                  |

## Conclusion

Monoisotopic precursor selection was demonstrated for TOF/TOF analysis of a standard TAG and TAGs in an olive oil sample. This selectivity made it possible to use charge-remote fragmentation to determine the complete structure (except stereochemistry) for all of the TAGs, including structural isomers, present in the sample. Multiple structural isomers in the precursor ion were identified through the observation of G- and J-type ions. These results demonstrate that the MALDI-TOF-TOF system with high precursor ion selectivity can fully analyze the structure of triacylglycerols without prior chromatographic separation, and that the method is effective for the analysis of complex fat composites in food.

## Open Access

This article is distributed under the terms of the Creative Commons Attribution License which permits any use, distribution, and reproduction in any medium, provided the original author(s) and the source are credited.

## References

- Byrdwell, W.C., Emken, E.A.: Analysis of triglycerides using atmospheric pressure chemical ionization mass spectrometry. *Lipids* **30**, 173–175 (1992)
- Neff, W.E., Byrdwell, W.C.: Soybean oil triacylglycerol analysis by reversed-phase high-performance liquid chromatography coupled with atmospheric pressure chemical ionization mass spectrometry. *J. Am. Oil Chem. Soc.* **72**, 1185–1191 (1995)
- Byrdwell, W.C., Emken, E.A., Neff, W.E., Adlof, R.O.: Quantitative analysis of triglycerides using atmospheric pressure chemical ionization mass spectrometry. *Lipids* **31**, 919–935 (1996)
- Byrdwell, W.C.: Atmospheric pressure chemical ionization mass spectrometry for the analysis of lipids. *Lipids* **36**, 327–346 (2001)
- Xu, Y., Brenna, J.T.: Atmospheric pressure covalent adduct chemical ionization tandem mass spectrometry for double bond location in monoene- and diene-containing triacylglycerols. *Anal. Chem.* **79**, 2525–2536 (2007)
- Hsu, F., Turk, J.: Structural characterization of triacylglycerols as lithiated adducts by electrospray ionization mass spectrometry using low-energy collisionally activated dissociation on a triple stage quadrupole instrument. *J. Am. Soc. Mass Spectrom.* **10**, 587–599 (1999)
- Byrdwell, W.C., Neff, W.E.: Dual parallel electrospray ionization and atmospheric pressure chemical ionization mass spectrometry (MS), MS/MS and MS/MS/MS for the analysis of triacylglycerols and triacylglycerol oxidation products. *Rapid Commun. Mass Spectrom.* **16**, 300–319 (2002)
- McAnoy, A., Wu, C., Murphy, R.: Direct qualitative analysis of triacylglycerols by electrospray mass spectrometry using a linear ion trap. *J. Am. Soc. Mass Spectrom.* **16**, 1498–1509 (2005)

9. Ham, B.M., Jacob, T.J., Keese, M.M., Cole, R.B.: Identification, quantification and comparison of major non-polar lipids in normal and dry eye tear lipidomes by electrospray tandem mass spectrometry. *J. Mass Spectrom.* **39**, 1321–1336 (2004)
10. Lévêque, N.L., Héron, S., Tchaplà, A.: Regioisomer characterization of triacylglycerols by non-aqueous reversed-phase liquid chromatography/electrospray ionization mass spectrometry using silver nitrate as a postcolumn reagent. *J. Mass Spectrom.* **45**, 284–296 (2010)
11. Gross, M.L.: Charge-remote fragmentations: method, mechanism and applications. *Int. J. Mass Spectrom. Ion Process.* **118/119**, 137–165 (1992)
12. Griffiths, W.J., Yang, Y., Lindgren, J.Å.: Charge Remote Fragmentation of Fatty Acid Anions in 400 eV Collisions with Xenon Atoms. *Rapid Commun. Mass Spectrom.* **10**, 21–28 (1996)
13. Deterding, J.D., Gross, M.L.: Tandem mass spectrometry for identifying fatty acid derivatives that undergo charge-remote fragmentations. *Org. Mass Spectrom.* **23**, 169–177 (1988)
14. Deterding, J.D., Gross, M.L.: Fast-atom-bombardment and tandem mass spectrometry for determining structures of fatty acids as their picolinyl ester derivatives. *Anal. Chim. Acta* **200**, 431–455 (1987)
15. Cordero, M.M., Wesdemiotis, C.: Characterization of the Neutral Products Formed upon Charge-Remote Fragmentation of Fatty Acid Ions. *Anal. Chem.* **66**, 861–866 (1994)
16. Murphy, R.C., Fiedler, J., Hevko, J.: Analysis of nonvolatile lipids by mass spectrometry. *Chem. Rev.* **101**(2), 479–526 (2001)
17. Ann, Q., Adams, J.: Structure determination of ceramides and neutral glycosphingolipids by collisional activation of  $[M + Li]^+$  ions. *J. Am. Soc. Mass Spectrom.* **3**, 260–263 (1992)
18. Trimpin, S., Clemmer, D.E., McEwen, C.N.: Charge-remote fragmentation of lithiated fatty acids on a TOF-TOF instrument using matrix-ionization. *J. Am. Soc. Mass Spectrom.* **18**, 1967–1972 (2007)
19. Jensen, N.J., Tomer, K.B., Gross, M.L.: FAB MS/MS for phosphatidylinositol, -glycerol, -ethanolamine, and other complex phospholipids. *Lipids* **22**, 480–489 (1987)
20. Shimma, S., Kubo, A., Satoh, T., Toyoda, M.: Detailed structural analysis of lipids directly on tissue specimens using a MALDI-SpiralTOF-Reflectron TOF mass spectrometer. *PLoS One* **7**, 5 (2012)
21. Cheng, C., Gross, M.L., Pittenauer, E.: Complete structural elucidation of triacylglycerols by tandem sector mass spectrometry. *Anal. Chem.* **70**, 4417–4426 (1998)
22. Cornish, T.J., Cotter, R.J.: Collision-induced dissociation in a tandem time-of-flight mass spectrometer with two single-stage reflectrons. *Org. Mass Spectrom.* **28**(10), 1129–1134 (1993)
23. Pittenauer, E., Allmaier, G.: The renaissance of high-energy CID for structural elucidation of complex lipids: MALDI-TOF/RTOF-MS of alkali cationized triacylglycerols. *J. Am. Soc. Mass Spectrom.* **20**, 1037–1047 (2009)
24. Satoh, T., Sato, T., Kubo, A., Tamura, J.: Tandem Time-of-Flight Mass Spectrometer with high precursor ion selectivity employing spiral ion trajectory and improved offset parabolic reflectron. *J. Am. Soc. Mass Spectrom.* **22**, 797–803 (2011)
25. Ishihara, M., Toyoda, M., Matsuo, T.: Perfect Spatial and Isochronous Focusing Ion Optics for Multi-turn Time of Flight Mass Spectrometer. *Int. J. Mass Spectrom.* **197**, 179–189 (2000)
26. Satoh, T., Tsuno, H., Iwanaga, M., Kammei, Y.: The design and characteristic features of a new time-of-flight mass spectrometer with a spiral ion trajectory. *J. Am. Soc. Mass Spectrom.* **16**(12), 1969–1975 (2005)

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
Fig.1 Structure of tristearin.


Fig.2 Mass spectrum of tristearin.


Fig.3 Product ion spectrum of sodium adducted tristearin.


Fig.4 Product ion spectrum of sodium adducted tristearin(enlarged between  $m/z$  650 and  $m/z$  920).


Fig.5 Peak assignment of obtained product ion spectrum.

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
Fig.1 Structure of triolein.


Fig.2 Mass spectrum of triolein.


Fig.3 Product ion spectrum of sodium adducted trierein.


Fig.4 Product ion spectrum of sodium adducted trierein(enlarged between  $m/z$  600 and  $m/z$  920).


Fig.5 Peak assignment of obtained product ion spectrum.

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
Fig.1 Structure of 1-Palmitoyl-2-oleoyl-3-linoleoyl-*rac*-glycerol and peak assignment of obtained product ion spectra.


Fig.2 Product ion spectrum of sodium adducted triacylglycerol.


Fig.3 Structure of J2 ion [3].


Fig.4 Product ion spectrum of sodium adducted triacylglycerol (enlarged between  $m/z$  600 and  $m/z$  920).

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
Figure.2. Mass spectrum of triolein after heating at 160C, 60min.

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
Figure 3. Product ion spectrum at  $m/z$  907 (top) and  $m/z$  923 (bottom)


Figure.4. Product ion spectrum at  $m/z$  923 (enlarged between  $m/z$  680 and  $m/z$  920).


Figure.5. Assignment of product ions at  $m/z$  923

[1] Cheng, C., Gross, M. L.; Pittenauer, E. Complete structural elucidation of triacylglycerols by tandem sector mass spectrometry. *Anal. Chem.* 1998, 70, 4417-4426

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
Fig.1 Chemical structure of phosphatidylcholine and peak assignment of obtained product ion spectra.


Fig.2 Mass spectrum of PC(16:0,18:1).


Fig.3 Product ion spectrum of PC(16:0,18:1) ( $[M+H]^+$ ).


Fig.4 Product ion spectrum of PC(16:0,18:1) (enlarged between  $m/z$  600 and  $m/z$  920).

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
Fig.1 Mass spectra of phospholipids from egg yolk (top:positive ion mode, bottom negative ion mode).


Fig.2 Product ion spectrum of PC(34:1) [M+H]<sup>+</sup>


Fig.3 Product ion spectrum of PI(38:4) [M-H]<sup>-</sup>


Fig.4 Peak assignment of product ion spectrum of PI(38:4) [M-H]<sup>-</sup>

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
Fig.1 Structure of 4ME 16:0 PC and peak assignment of obtained product ion spectra.


Fig.2 Mass spectrum of 4ME 16:0 PC.


Fig.3 Product ion spectrum of 4ME 16:0 PC  $[M+H]^+$ .


Fig.4 Product ion spectrum of 4ME 16:0 PC  $[M+H]^+$  (enlarged between  $m/z$  500 and  $m/z$  850).

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
Fig.1 Structure of GM1 and peak assignment of obtained product ion spectra.


Fig.2 Mass spectrum of GM1.


Fig.3 Product ion spectrum of  $[M_1-H]^-$ .


Fig.4 Product ion spectrum of  $[M_1-H]^-$  (enlarged between  $m/z$  1120 and  $m/z$  1540).


Fig.5 Structure of Ceramide moiety of GM1 and peak assignment of obtained product ion spectra.

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
Fig.1 Product ion spectrum of renin substrate tetradecapeptide.


Fig.2 Product ion spectrum of renin substrate tetradecapeptide and Pathway for d ion (Left: Isoleucine, Right: Leucine).


Fig.3 Structure of immonium ion (R: amino-acid residue)


Fig.4 Product ion spectrum of renin substrate tetradecapeptide enlarged between  $m/z$  50 and 150.

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
Fig.1 Mass spectrum of Monophosphopeptide.


Fig.2 Product ion spectrum of Monophosphopeptide.

- [1] A. Stensballe, O.N. Jensen, J.V. Olsen, K.F. Haselmann, R.A. Zubarev, Rapid Communications in Mass Spectrometry 14 (2000) 1793.
- [2] S. Shimma, H. Nagao, A.E. Giannakopoulos, S. Hayakawa, K. Awazu, M. Toyoda, Journal of Mass Spectrometry 43 (2008) 535.

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Table1 Experimental conditions and specifications

|                   | JMS-S3000     | JMS-700T      |
|-------------------|---------------|---------------|
| Ionization method | MALDI         | FAB (Xe 6 kV) |
| Sample quantity   | 2.5 pmol      | 1 nmol        |
| Precursor ion     | mono isotopic | mono isotopic |
| Target gas        | He            | Ar            |
| Matrix            | CHCA          | Glycerin, NBA |
| Collision energy  | 20000 eV      | 5000 eV       |


Fig.1 Ion trajectory Left: "JMS-S3000 SpiralfTOF", Right: "JMS-700T tandem MStation"


Fig.2 Product ion spectra of Renin-substrate tetradecapeptide by  
 a) JMS-S3000 SpiralfTOF-TOF and b) JMS-700T tandem MStation

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
Fig.1 Mass spectrum of SAAF.


Fig.2 Product ion spectrum of  $[M-2Na-SO_3+H]^+$


Fig.3 Structure of SAAF and fragmentation pattern (Asterisks show dehydrated ions.)

[1] M. Yoshida, M. Murata, K. Inaba, and M. Morisawa, Proceeding of the National Academy of Sciences of the United States of America 99 (2002) 14831

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
Fig.1 Mass spectrum of YTX.


Fig.2 Product ion spectrum of  $[M-2Na-SO_3+H]^+$


Fig.3 Structure of YTX and fragmentation pattern

[1] H. Naoki, M. Murata, T. Yasumoto, Rapid Communication of Mass Spectrometry 7 (1993) 179

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- 32 Sato, H., Ishii, Y., Momose, H., Sato, T. & Teramoto, K. Application of High-Resolution MALDI-TOFMS with a Spiral Ion Trajectory for the Structural Characterization of Free Radical Polymerized Methacrylate Ester Copolymers. *Mass Spectrom (Tokyo)* **2**, A0014, doi:10.5702/massspectrometry.A0014 (2013).
- 33 Teramoto, K. *et al.* Simple and rapid characterization of mycolic acids from Dietzia strains by using MALDI spiral-TOFMS with ultra high mass-resolving power. *The Journal of antibiotics* **66**, 713-717, doi:10.1038/ja.2013.79 (2013).
- 34 Tsujita, T., Shintani, T. & Sato, H.  $\alpha$ -Amylase inhibitory activity from nut seed skin polyphenols. 1. Purification and characterization of almond seed skin polyphenols. *Journal of agricultural and food chemistry* **61**, 4570-4576, doi:10.1021/jf400691q (2013).
- 35 Voorhees, K. J. *et al.* Modified MALDI MS fatty acid profiling for bacterial identification. *Journal of Mass Spectrometry* **48**, 850-855, doi:10.1002/jms.3215 (2013).
- 36 Wang, X., Fukuoka, S., Tsukigawara, R., Nagata, K. & Higuchi, M. Electric-field-enhanced oriented cobalt coordinated peptide monolayer and its electrochemical properties. *Journal of colloid and interface science* **390**, 54-61, doi:10.1016/j.jcis.2012.08.079 (2013).
- 37 Yamada, S., Takasu, A. & Kawamura, K. The effect of microwave irradiation on the kinetics and activation thermodynamics of ring - opening polymerization of  $\epsilon$  - caprolactone. *Journal of Polymer Science Part A: Polymer Chemistry* **51**, 3732-3739, doi:10.1002/pola.26776 (2013).
- 38 Akabane-Nakata, M., Obika, S. & Hari, Y. Synthesis of oligonucleotides containing N, N-disubstituted 3-deazacytosine nucleobases by post-elongation modification and their triplex-forming ability with double-stranded DNA. *Organic & biomolecular chemistry* **12**, 9011-9015, doi:10.1039/C4OB01760A (2014).
- 39 Chen, Y.-T., Hong, P.-F., Wen, L. & Lin, C.-T. Molecular cloning and characterization of a thioredoxin from Taiwanofungus camphorata. *Botanical Studies* **55**, 1-8, doi:10.1186/s40529-014-0077-z (2014).
- 40 Ito, M., Shimasaki, T., Teramoto, N. & Shibata, M. Photo-cured organic-inorganic hybrid composites of methacrylate-terminated 4-arm star-shaped  $\epsilon$ -caprolactone oligomer and methacrylate-or thiol-substituted polysilsesquioxane. *Journal of Polymer Research* **21**, 1-10, doi:10.1007/s10965-014-0507-3 (2014).
- 41 Kawasaki, T., Imai, T. & Tsukiyama, K. Use of a Mid-Infrared Free-Electron Laser (MIR-FEL) for Dissociation of the Amyloid Fibril Aggregates of a Peptide. *Journal of Analytical Sciences, Methods and Instrumentation* **04**, 9-18, doi:10.4236/jasmi.2014.41002 (2014).
- 42 Ken, C.-F., Chang, C.-C., Wen, L., Huang, J.-K. & Lin, C.-T. Biochemical characterization of a functional recombinant aryl-alcohol dehydrogenase from Taiwanofungus camphorata. *Botanical Studies* **55**, 14, doi:10.1186/1999-3110-55-14 (2014).
- 43 Kim, B.-J. *et al.* Preparation of carbon fibers with excellent mechanical properties from isotropic pitches. *Carbon* **77**, 747-755, doi:10.1016/j.carbon.2014.05.079 (2014).
- 44 Kurashige, W. *et al.* Au<sub>25</sub> Clusters Containing Unoxidized Tellurolates in the Ligand Shell. *The journal of physical chemistry letters* **5**, 2072-2076, doi:10.1021/jz500901f (2014).
- 45 Niihori, Y., Matsuzaki, M., Uchida, C. & Negishi, Y. Advanced use of high-performance liquid chromatography for synthesis of controlled metal clusters. *Nanoscale* **6**, 7889-7896, doi:10.1039/C4NR01144A (2014).
- 46 Nishikawa, K. *et al.* Detection of changes in the structure and distribution map of triacylglycerol in fatty liver model by MALDI - SpiralTOF. *FEBS open bio* **4**, 179-184, doi:10.1016/j.fob.2014.02.005 (2014).
- 47 Ohtani, H. & Iura, T. Complementary Characterization of End Groups in Radically Polymerized Poly(methyl methacrylate) by Matrix-Assisted Laser Desorption/Ionization Mass Spectrometry and Pyrolysis-Gas Chromatography-Mass Spectrometry. *Mass Spectrom (Tokyo)*

- 3, S0041, doi:10.5702/massspectrometry.S0041 (2014).
- 48 Okabe, T., Kuzuhara, D., Aratani, N. & Yamada, H. Synthesis and electronic properties of acetylene-and butadiyne-linked 3, 3'-porphycene dimers. *Journal of Porphyrins and Phthalocyanines* **18**, 849-855, doi:10.1142/S1088424614500631 (2014).
- 49 Sanai, Y., Morita, Y., Asano, Y., Ishizaki, K. & Kubota, K. Chain - end lactonization of polyacrylates prepared by photopolymerization. *Journal of Polymer Science Part A: Polymer Chemistry* **52**, 1161-1171, doi:10.1002/pola.27100 (2014).
- 50 Sato, H., Nakamura, S., Teramoto, K. & Sato, T. Structural characterization of polymers by MALDI spiral-TOF mass spectrometry combined with Kendrick mass defect analysis. *J Am Soc Mass Spectrom* **25**, 1346-1355, doi:10.1007/s13361-014-0915-y (2014).
- 51 Satoh, T., Kubo, A., Hazama, H., Awazu, K. & Toyoda, M. Separation of Isobaric Compounds Using a Spiral Orbit Type Time-of-Flight Mass Spectrometer, MALDI-SpiralTOF. *Mass Spectrom (Tokyo)* **3**, S0027, doi:10.5702/massspectrometry.S0027 (2014).
- 52 Shibata, T., Fujimoto, M. & Otani, T. Synthesis of macrocyclic heteroarylenes by consecutive inter-and intramolecular cycloadditions of thiophenylene-tethered triynes. *Tetrahedron* **70**, 8453-8461, doi:10.1016/j.tet.2014.09.077 (2014).
- 53 Shigeri, Y. *et al.* Hydrazide and hydrazine reagents as reactive matrices for MALDI - MS to detect gaseous aldehydes. *Journal of Mass Spectrometry* **49**, 742-749, doi:10.1002/jms.3408 (2014).
- 54 Suzuki, R. *et al.* Brevisulcatic acids, marine ladder-frame polyethers from the red tide dinoflagellate *Karenia brevisulcata* in New Zealand. *Org Lett* **16**, 5850-5853, doi:10.1021/ol502700h (2014).
- 55 Takagi, K., Kawagita, E. & Kouchi, R. Synthesis and characterization of polythiophene derivatives with nitrogen heterocycles on the side chain. *Journal of Polymer Science Part A: Polymer Chemistry* **52**, 2166-2174, doi:10.1002/pola.27228 (2014).
- 56 Tsujita, T., Shintani, T. & Sato, H. Preparation and characterisation of peanut seed skin polyphenols. *Food chemistry* **151**, 15-20, doi:10.1016/j.foodchem.2013.11.072 (2014).
- 57 Wambua, D. M., Ubukata, M., Dane, J., Cody, R. B. & Chiu, N. H. Bottom-up mass spectrometric sequencing of microRNA. *Analytical Methods* **6**, 8829-8839, doi:10.1039/C4AY01519C (2014).
- 58 Yamaguchi, I. & Nagano, T. Synthesis, chemical, and thermoelectric properties of n - type  $\pi$  - conjugated polymer composed of 1, 2, 4 - triazole and pyridine rings and its metal complexes. *Journal of Applied Polymer Science* **131**, doi:10.1002/app.39928 (2014).
- 59 Yamakado, R. *et al.* Synthesis, reaction, and optical properties of cyclic oligomers bearing 9, 10-diphenylanthracene based on an aromatic tertiary amide unit. *RSC Advances* **4**, 6752-6760, doi:10.1039/C3RA46652C (2014).
- 60 Yamazoe, S., Kurashige, W., Nobusada, K., Negishi, Y. & Tsukuda, T. Preferential Location of Coinage Metal Dopants (M= Ag or Cu) in  $[\text{Au}_{25-x}\text{M}_x(\text{SC}_2\text{H}_4\text{Ph})_{18}]^{-(x-1)}$  As Determined by Extended X-ray Absorption Fine Structure and Density Functional Theory Calculations. *The Journal of Physical Chemistry C* **118**, 25284-25290, doi:10.1021/jp5085372 (2014).
- 61 Yan, Y., Ubukata, M., Cody, R. B., Holy, T. E. & Gross, M. L. High-energy collision-induced dissociation by MALDI TOF/TOF causes charge-remote fragmentation of steroid sulfates. *J Am Soc Mass Spectrom* **25**, 1404-1411, doi:10.1007/s13361-014-0901-4 (2014).
- 62 Yu, M., Tang, T., Takasu, A. & Higuchi, M. pH-and thermo-induced morphological changes of an amphiphilic peptide-grafted copolymer in solution. *Polymer journal* **46**, 52-58, doi:10.1038/pj.2013.61 (2014).
- 63 Abe, T., Matsuzaka, M., Isoda, K. & Tadokoro, M. Liquid-Crystalline and Electronic Properties

- of Racemic-Alkoxy Chains-Substituted Tetraazanaphthacene. *Molecular Crystals and Liquid Crystals* **615**, 70-77, doi:10.1080/15421406.2015.1066963 (2015).
- 64 Asamizu, S., Ozaki, T., Teramoto, K., Satoh, K. & Onaka, H. Killing of Mycolic Acid-Containing Bacteria Aborted Induction of Antibiotic Production by Streptomyces in Combined-Culture. *PLoS One* **10**, e0142372, doi:10.1371/journal.pone.0142372 (2015).
- 65 Hayakawa, T., Song, H., Ishii, Y. & Kawasaki, S. Facile Synthesis of Coronene Dimers Having Unique Optical Properties. *J. Chem* **9**, 442-447, doi:10.17265/1934-7375/2015.07.002 (2015).
- 66 Hayakawa, Y. *et al.* Complementary DNA display selection of high - affinity peptides binding the vacuolating toxin (VacA) of Helicobacter pylori. *Journal of Peptide Science* **21**, 710-716, doi:10.1002/psc.2795 (2015).
- 67 Hubin, E. A. *et al.* Structural, functional, and genetic analyses of the actinobacterial transcription factor RbpA. *Proc Natl Acad Sci U S A* **112**, 7171-7176, doi:10.1073/pnas.1504942112 (2015).
- 68 Ikawa, T. *et al.* 2 - [(Neopentyl glycolato) boryl] phenyl Triflates and Halides for Fluoride Ion - Mediated Generation of Functionalized Benzynes. *Advanced Synthesis & Catalysis* **357**, 2287-2300, doi:10.1002/adsc.201500315 (2015).
- 69 Itoh, Y. *et al.* Identification of SNAIL1 Peptide-Based Irreversible Lysine Specific Demethylase 1-Selective Inactivators. *Journal of medicinal chemistry*, doi:10.1021/acs.jmedchem.5b01323 (2015).
- 70 Iura, T. & Ohtani, H. Fragmentation behavior of poly(methyl methacrylate) during matrix-assisted laser desorption/ionization. *Rapid Commun Mass Spectrom* **29**, 155-162, doi:10.1002/rm.7087 (2015).
- 71 Kano, H. *et al.* A novel analytical pyrolysis device applicable for measurements of less volatile pyrolyzates. *Journal of Analytical and Applied Pyrolysis* **113**, 165-173, doi:10.1016/j.jaap.2014.12.008 (2015).
- 72 Kurimoto, K., Yamamura, H. & Miyagawa, A. Chemical approach for the syntheses of GM4 isomers with sialic acid to non-natural linkage positions on galactose. *Carbohydrate research* **401**, 39-50, doi:10.1016/j.carres.2014.10.018 (2015).
- 73 Miyagawa, A., Matsuda, T. & Yamamura, H. Synthesis of Branched Tetrasaccharide Derivatives of Schizophyllan-like  $\beta$ -Glucan. *Journal of Carbohydrate Chemistry* **34**, 215-246, doi:10.1080/07328303.2015.1044754 (2015).
- 74 Mizukado, J. *et al.* High - resolution MALDI - TOF MS study on analysis of low - molecular - weight products from photo - oxidation of poly (3 - hexylthiophene). *Journal of Mass Spectrometry* **50**, 1006-1012, doi:10.1002/jms.3614 (2015).
- 75 Mori, S., Morihiro, K., Kasahara, Y., Tsunoda, S.-i. & Obika, S. Synthesis and Properties of 2'-Deoxyuridine Analogues Bearing Various Azobenzene Derivatives at the C5 Position. *Chemosensors* **3**, 36-54, doi:10.3390/chemosensors3020036 (2015).
- 76 Negishi, Y. *et al.* A critical size for emergence of nonbulk electronic and geometric structures in dodecanethiolate-protected Au clusters. *Journal of the American Chemical Society* **137**, 1206-1212, doi:10.1021/ja5109968 (2015).
- 77 Niihori, Y., Kikuchi, Y., Kato, A., Matsuzaki, M. & Negishi, Y. Understanding Ligand-Exchange Reactions on Thiolate-Protected Gold Clusters by Probing Isomer Distributions Using Reversed-Phase High-Performance Liquid Chromatography. *ACS nano* **9**, 9347-9356, doi:10.1021/acsnano.5b03435 (2015).
- 78 Nishikawa, K. *et al.* Resveratrol increases CD68+ Kupffer cells colocalized with adipose differentiation - related protein and ameliorates high - fat - diet - induced fatty liver in mice. *Molecular nutrition & food research* **59**, 1155-1170, doi:10.1002/mnfr.201400564 (2015).

- 79 Ono, K. *et al.* Self-Assembly of Nanometer-Sized Boroxine Cages from Diboronic Acids. *Journal of the American Chemical Society* **137**, 7015-7018, doi:10.1021/jacs.5b02716 (2015).
- 80 Sakano, T., Ohashi, T., Yamanaka, M. & Kobayashi, K. Photoresponsive self-assembled hexameric capsules based on calix [4] resorcinarenes bearing azobenzene dendron conjugates as side chains. *Organic & biomolecular chemistry* **13**, 8359-8364, doi:10.1039/C5OB00997A (2015).
- 81 Satoh, T. *et al.* Imaging Mass Spectrometry Using Ultra-high Mass Resolution Matrix-Assisted Laser Desorption/Ionization Time-of-Flight Mass Spectrometer, SpiralTOF. *Microscopy and Microanalysis* **21**, 2059-2060, doi:10.1017/S1431927615011071 (2015).
- 82 Segawa, Y., Kuwabara, T., Matsui, K., Kawai, S. & Itami, K. Palladium-free synthesis of [10] cycloparaphenylene. *Tetrahedron* **71**, 4500-4503, doi:10.1016/j.tet.2015.02.066 (2015).
- 83 Sharma, S., Kurashige, W., Nobusada, K. & Negishi, Y. Effect of trimetallization in thiolate-protected Au(24-n)Cu(n)Pd clusters. *Nanoscale* **7**, 10606-10612, doi:10.1039/c5nr01491c (2015).
- 84 Shigeri, Y. *et al.* Hydrazide and hydrazine reagents as reactive matrices for matrix-assisted laser desorption/ionization mass spectrometry to detect steroids with carbonyl groups. *European Journal of Mass Spectrometry* **21**, 79-90, doi:10.1255/ejms.1336 (2015).
- 85 Spring, A. M. *et al.* An analysis of the structural, thermal and optical characteristics as well as the electrical resistivity of tert-butyldiphenylsilyl substituted poly (norbornene-dicarboximide)s. *Polymer* **56**, 189-198, doi:10.1016/j.polymer.2014.11.043 (2015).
- 86 Takanashi, E. *et al.* Cytotoxic linear acetylenes from a marine sponge Pleroma sp. *Tetrahedron* **71**, 9564-9570, doi:10.1016/j.tet.2015.10.062 (2015).
- 87 Teramoto, K. *et al.* Characterization of Mycolic Acids in Total Fatty Acid Methyl Ester Fractions from Mycobacterium Species by High Resolution MALDI-TOFMS. *Mass Spectrom (Tokyo)* **4**, A0035, doi:10.5702/massspectrometry.A0035 (2015).
- 88 Wang, C. *et al.* A novel D-π-A small molecule with N-heteroacene as acceptor moiety for photovoltaic application. *Dyes and Pigments* **122**, 231-237, doi:10.1016/j.dyepig.2015.06.029 (2015).
- 89 Wang, C. *et al.* Synthesis, Characterization, and Memory Performance of Two Phenazine/Triphenylamine - Based Organic Small Molecules through Donor - Acceptor Design. *Asian Journal of Organic Chemistry* **4**, 646-651, doi:10.1002/ajoc.201500087 (2015).
- 90 Yamanoi, T., Oda, Y., Katsuraya, K., Inazu, T. & Hattori, K. Synthesis, structure, and evaluation of a β-cyclodextrin-artificial carbohydrate conjugate for use as a doxorubicin-carrying molecule. *Bioorganic & medicinal chemistry*, doi:10.1016/j.bmc.2015.12.030 (2015).
- 91 Zheng, Q., Morimoto, M., Sato, H. & Takanohashi, T. Molecular composition of extracts obtained by hydrothermal extraction of brown coal. *Fuel* **159**, 751-758, doi:10.1016/j.fuel.2015.07.042 (2015).
- 92 橋詰峰雄 *et al.* 体液類似環境におけるアパタイト結合性ペプチド樹状オリゴマーによるミネラルイオン. *高分子論文集* **73**, 55-61, doi:10.1295/koron.2015-0050 (2015).
- 93 小林大祐 *et al.* 超音波を用いたポリピロール微粒子合成. *化学工学論文集* **41**, 153-156, doi:10.1252/kakoronbunshu.41.153 (2015).
- 94 Fouquet, T., Nakamura, S. & Sato, H. MALDI SpiralTOF high-resolution mass spectrometry and Kendrick mass defect analysis applied to the characterization of poly(ethylene-co-vinyl acetate) copolymers. *Rapid Commun Mass Spectrom* **30**, 973-981, doi:10.1002/rcm.7525 (2016).
- 95 Hayakawa, T., Ishii, Y. & Kawasaki, S. Sodium ion battery anode properties of designed graphene-layers synthesized from polycyclic aromatic hydrocarbons. *RSC Advances* **6**,

22069-22073, doi:10.1039/C6RA00955G (2016).

- 96 Kakizawa, T. *et al.* Evaluation of phenylcyclopropylamine compounds by enzymatic assay of lysine-specific demethylase 2 in the presence of NPAC peptide. *Bioorganic & Medicinal Chemistry Letters*, doi:10.1016/j.bmcl.2016.01.036 (2016).
- 97 Katagiri, K., Takasu, A. & Higuchi, M. Synthesis of Glycopolymer Containing Cell-Penetrating Peptide as Inducers of Recombinant Protein Expression under the Control of Lac Operator/Repressor Systems. *Biomacromolecules*, doi:10.1021/acs.biomac.6b00368 (2016).
- 98 Kim, B.-J. *et al.* Enhancing the tensile strength of isotropic pitch-based carbon fibers by improving the stabilization and carbonization properties of precursor pitch. *Carbon* **99**, 649-657, doi:10.1016/j.carbon.2015.12.082 (2016).
- 99 Peršurić, Ž., Osuga, J., Grbac, T. G., Peter - Katalinić, J. & Pavelić, S. K. MALDI - SpiralTOF technology for assessment of triacylglycerols in Croatian olive oils. *European Journal of Lipid Science and Technology*, doi:10.1002/ejlt.201500375 (2016).
- 100 Shigeri, Y. *et al.* 2-hydrazinoquinoline: a reactive matrix for MALDI-MS to detect gaseous carbonyl compounds. *European Journal of Mass Spectrometry* **22**, 0-0, doi:10.1255/ejms.1413 (2016).
- 101 Takahira, Y., Murotani, E., Fukuda, K., Vohra, V. & Murata, H. Design, synthesis, and properties of a series of charged iridium (III) complexes with a neutral bidentate ligand for deep-blue phosphorescent emitter. *Journal of Fluorine Chemistry* **181**, 56-60, doi:10.1016/j.jfluchem.2015.10.016 (2016).
- 102 Yang, J., Nakabayashi, K., Miyawaki, J. & Yoon, S.-H. Preparation of isotropic pitch-based carbon fiber using hyper coal through co-carbonation with ethylene bottom oil. *Journal of Industrial and Engineering Chemistry* **34**, 397-404, doi:10.1016/j.jiec.2015.11.026 (2016).

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Copyright(C)2016 JEOL Ltd. All Rights Reserved.

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