

Utilities of No-D NMR (No-Deuterium Proton NMR)

Product used : Nuclear Magnetic Resonance (NMR)

No-D NMR (No-Deuterium Proton NMR) technique is a measurement of high resolution ¹H NMR spectra without deuterium solvent. It suggests that any reaction mixture or reagent solution are directly available for No-D NMR measurement.

In conventional NMR measurements, ²H signals of deuterium solvent are used for the shimming, however ¹H solvent signals are as well in No-D NMR measurements. Since suppression of strong ¹H solvent signals by WET eliminates ¹³C satellite signals also, it is convenient approach to collect ¹H NMR spectra without deuterium solvent by No-D NMR measurement.

Merits of No-D NMR

“No-D NMR “is your solution!

- The better approach for unstable compound just after chemical reactions
- Skipping the laborious sample preparation (dissolving into a deuterated solvent after sample purification)
- Cost reduction for deuterated solvent in daily NMR measurement work

Automation script for No-D NMR

The following 5 steps are automated No-D NMR measurement procedure on NMR software Delta, and every tedious measurement parameters are automatically adjusted.

- ① Shimming
- ② Detecting solvent signals
- ③ Suppression of solvent signals
- ④ Processing FID data
- ⑤ Chemical shift adjustment

On Delta software, chemical shift in the conventional ¹H measurement is not adjusted without NMR lock (Fig. 1a), but the automation of No-D NMR adjusts chemical shift automatically (Fig. 1b). Its result has the same chemical shift information as deuterated solvent sample spectrum (Fig. 1c).




Fig.1 ¹H spectra of 5 mg rotenone measured with JNM-ECZ500R

For starting No-D NMR measurement, only simple 2 parameters are required (Fig.2).

1. Select solvent
2. Input the number of suppression signals




Fig.2 Operation window of No-D NMR (Delta V5)

Application of No-D NMR 1: Observation of exchangeable ¹H signals

When the compound with exchangeable proton is dissolved in exchangeable proton solvent (D₂O, CD₃OD, and etc.), proton atom is replaced with deuterium atom. As results, it is not possible to observe amino protons for the target molecule (Fig. 3 upper). On the other hand, the exchangeable proton can be observed in No-D NMR, since amino protons are kept under protonated solvent (Fig. 3 bottom).




Fig.3 2.5mg L-Tryptophan measured with JNM-ECZ600R

Application of No-D NMR 2: Apply to mixed solvent

Because No-D NMR can suppress multi site signals, it is also effective for multi signal solvents such as CH₃OH and mixed solvent sample. Furthermore, ¹³C decoupling is useful, when solvent ¹³C satellite signal overlap with sample signal. Fig.5 indicates how sample signal can be distinguished from ¹³C satellite by ¹³C decoupling. Chemical shifts in Fig.4 and Fig.5 are intentionally adjusted for better clarity.




Fig.5 expansion of Fig.4

Fig.4 2.5mg ibuprofen (racemic form)
Measured with JNM-ECZ400S

