KNT. R2111HD/NI Application Data

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Table 1 HS/Sniffing GC/MS measurement conditions.

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
Fig.1 TIC chromatogram of Earl Grey 5 minute extraction sample.


Fig.2 TIC chromatogram of Earl Grey leaf sample.

Table 2 Retention time, identified compounds and smell in Earl Grey 5 minute extraction sample.

Number	R.T.	Estimated compound name	Smell	Intensity
1	2:54	Dimethylsulfide	Aroma	Middle
2	3:44	2-Methylbutanal	Grass	Middle
3	4:49	Pinene	Bond	Weak
4	5:49	Nopinene	Green	Middle
5	7:01	Limonene	Mint	Strong
6	7:40	cis-Ocimene	Grass	Weak
7	7:59	2-Carene	Flower	Middle
8	9:13	2-methyl-2-hepten-6-one	Flower	Middle
9	11:01	Furfural	Cheese	Middle
10	12:15	Linalool	Rose	Strong
11	13:08	Phenylacetaldehyde	Green	Middle
12	13:38	(Z)-2-(3,3-Dimethylcyclohexylidene)ethanol	Aroma	Middle
13	13:51	α -Terpineol	Flower	Middle
14	15:25	Geraniol	Flavor	Middle

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Table 3 Retention time, identified compounds and smell in Earl Grey leaf sample.

Number	R.T.	Estimated compound name	Smell	Intensity
1	2:56	Dimethylsulfide	Flavor	Strong
2	3:51	2-Methylbutanal	Grass	Middle
3	5:14	Pinene	Bond	Middle
4	6:14	Nopinene	Grass	Strong
5	7:36	Limonene	Mint	Strong
6	8:18	Terpinene	Grassy mint	Middle
7	8:27	p-Cymene	Grassy mint	Middle
8	8:37	2-Carene	Flower	Middle
9	9:20	2-Methyl-2-hepten-6-one	Fruit	Middle
10	9:44	Nealloocimene	Flavor	Middle
11	9:50	2,4-Dimethylpyrazine	Flavor	Middle
12	10:09	cis-3-Hexanol	Fresh green	Middle
13	11:01	Furfural	Baked cheese	Middle
14	12:16	Linalool	Rose	Strong
15	12:31	Linalylanthranilate	Fresh orange	Strong
16	13:11	Phenylacetaldehyde	Green	Middle
17	13:51	α -Terpineol	Flower	Middle
18	15:25	Geraniol	Sweet grass	Middle
19	16:05	Hydroxyethylbenzene	Delicate mint	Weak
20	16:41	2-Acetylpyrrole	Sweet	Middle
21	17:16	Pantolactone	Sweet	Strong
22	19:28	3,5-Dihydroxy-6-methyl-2H-pyran-4(3H)-one	Grape	Strong


Fig.3 Correlation between a extraction time and the peak intensity of smell compounds in Earl Grey sample.