



Solid-State NMR Meets Electron Diffraction

To distinguish crystalline form of small molecule, powder X-ray diffraction and ^{13}C CPMAS solid-state NMR are widely used. Here, we introduce a combined approach of ^1H solid-state NMR at ultrafast MAS and electron diffraction using TEM to overcome the difficulty which powder X-ray and ^{13}C CPMAS solid-state NMR have.



JEOL's philosophy for open innovation!
Here we combine solid-state NMR and TEM.



Copyright © 2017 JEOL Ltd.

Reference
Oikawa, T., Okumura, M., Kimura, T. & Nishiyama, Acta Cryst. C73 (2017) 219–228.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.